




# **RANSCREENER<sup>®</sup>** Assays

## Targeting the cGAS-STING Pathway Using a Homogenous, HTS Compatible Transcreeper<sup>®</sup> cGAS Assay

Meera Kumar, Justin Brink, and Robert G. Lowery  
BellBrook Labs, Madison, WI

### Quick Read

Poster presented September 2020 at  
SLAS in San Diego, CA




### Overview

Cyclic GMP-AMP synthase (cGAS) is a recently discovered enzyme that acts as a foreign DNA sensor to elicit an immune response to pathogens via activation of the STING (stimulator of interferon genes) receptor. Shortly after its discovery in 2013, aberrant activation of cGAS by self-DNA was shown to underlie debilitating and sometimes fatal autoimmune diseases, such as systemic lupus erythematosus (SLE) and Aicardi-Goutieres Syndrome (AGS). Knockout studies in animal models have clearly indicated that inhibiting cGAS is a promising approach for therapeutic intervention. However, there are no HTS-compatible assay methods for measuring cGAS enzyme activity or for monitoring cGAMP in cell and tissue samples. To enable HTS efforts targeting cGAS, we developed Transcreeper-based assays for cGAMP detection; i.e., homogenous, competitive immunoassays with fluorescence polarization (FP) and time resolved Förster resonance energy transfer (TR-FRET) signals. The key assay reagents are antibodies that selectively bind cGAMP in the presence of excess ATP and GTP and tracers that are displaced by cGAMP. We validated both assays with full length human cGAS and performed pilot screens for cGAS inhibitors with the FP assay. Similar to other Transcreeper assays, the cGAMP assays have the performance characteristics needed for small molecule screening including sensitivity, robustness (good dynamic range, minimal signal variation), and low levels of compound interference. The availability of robust, HTS compatible assays for measuring cGAS enzyme activity will accelerate efforts to target the cGAS-STING pathway for autoimmune diseases.




# The cGAS-cGAMP-STING Pathway Activates the Immune System in Response to Cytosolic DNA




**Figure 1.** IFN-driven immune responses triggered by cGAS/STING are critical for protection against microbial pathogens and for tumor cell-specific T cell responses in cancer, but activation by self-nucleic acids can contribute to autoimmune diseases.

## Mix and Read HTS Enzymatic Assays for cGAS Based on Immunodetection of cGAMP



**Figure 2.** Transcreener cGAS Assay principle: Enzymatically generated cGAMP displaces a fluorescent tracer from mAb causing a decrease in its polarization or TR-FRET.


## Assay Validation



**Figure 3. A. Specificity of mAb:** Measuring cGAS initial activity requires detection of cGAMP in an excess of ATP and GTP. Competition binding curves show outstanding selectivity for cGAMP vs. substrates, ATP and GTP, as well as related molecules.


**B. Detection of purified, full-length human cGAS with FP Assay.** cGAS enzyme reactions contained 100  $\mu$ M ATP and GTP, 62.5 nM 45 bp ISD DNA, 60 min reactions. N- and C-terminal His-tagged cGAS was produced at BBL; 6xHis-cGAS was also generously supplied by Z. Chen (UTSW Medical Center). C. Linear response: Polarization data from A. was converted to cGAMP using a standard curve. **D. Detection of purified, full-length human cGAS with TR-FRET Assay.** TR-FRET was used as an orthogonal assay method for hit confirmation. As with the FP assay, conversion of FRET signal to cGAMP formation yielded a linear response (data not shown). **E. Interference screen:** The LOPAC library (1280 compounds) was used at 10  $\mu$ M final concentration. Wells contained all cGAS enzyme reaction components except the cGAS enzyme; positive controls contained 1  $\mu$ M cGAMP, negative controls contained no cGAMP. **F. Pilot screen with 1600 diversity compounds (Life Chemicals):** cGAS was used at 10 nM, compounds were at 10  $\mu$ M; 60 min reaction; negative controls lacked dsDNA (required for cGAS activation);  $Z = 0.62$ ,  $Z' = 0.7$ .

## Screening and Triaging of Non-Stoichiometric Inhibitors




**Figure 4. A. HTS Workflow:** Primary screen and follow up assays used to triage undesirable compounds and select cGAS inhibitors for advancement into medicinal chemistry/SAR. **B. Scatter plot from 3,200 compounds (10 plates) from a 100K screen:** compounds at 20  $\mu$ M; negative controls lacked dsDNA (required for cGAS activation);  $Z = 0.59$ ,  $Z' = 0.63$ . **C., D. Confirmation of a screening hit using FP assay (C) and TR-FRET assay (D) for dose response.** Average IC<sub>50</sub>s were 12.8  $\mu$ M and 22.7  $\mu$ M, respectively. **E., F., G. H. Identification of non-stoichiometric inhibitors (NSI) using inhibitor titration (E., F.) and detergent disaggregation (G. H.).** Excess enzyme titrates out NSIs resulting in a significant increase in IC<sub>50</sub> (>3-fold), (F.). Non-ionic detergent can disperse aggregates, resulting in decreased IC<sub>50</sub> (H.).

## Characterizing Mode of Inhibition



**Figure 5. Characterization of Inhibitor MOA During Hit-to-Lead Using Substrate Competition Assays. A., B. Identification of a GTP-competitive inhibitor. C., D. Identification of an inhibitor that competes with both ATP and GTP. The Transcreeper cGAS FP assay was used for all assays.**

## Measuring Inhibitor Residence Time



**Figure 6. A. Schematic of jump dilution method used to measure inhibitor residence time. B. Monitoring recovery of activity for a series of cell-permeable inhibitors following jump dilution. Activity recovers as the E-I complex dissociates, allowing calculation of off-rates.**

## Conclusions

- Competitive immunoassays for cGAMP with FP and TR-FRET readouts were validated for detection of cGAS enzyme activity and small molecule screening.
- The assays served as a critical tool in a cGAS HTS campaign, including screening, hit confirmation and triaging of non-stoichiometric inhibitors.
- The Transcreeper cGAS FP assay was used to sort analogs based on mechanism of inhibition during hit-to-lead, which provided valuable information for driving SAR.
- The ability to run the assay in continuous mode allowed measurement of inhibitor residence time in an automated format, allowing incorporation of binding kinetics into SAR.