


INFINITI

Duality Vet Stent®



Study:

Evaluating Duality Vet Stent® in Treating Tracheal Collapse

- Retrospective review of 20 cases of tracheal collapse treated with Duality Vet Stent®
- Follow-up period of up to 17 months
- Zero stent fractures reported and only 1 out of 20 cases reported any granulation tissue
- Statistically significant improvement in cough, quality of life, exercise tolerance, and breathing
- 95% of clients would repeat the procedure and 90% were satisfied with the procedure

Score	Pre-Stent	Post-Stent	P-value
Cough 0 (none) to 10 (worst)	6.55 ± 2.44	2.5 ± 1.54	< 0.0001
Honking 0 (none) to 10 (worst)	6 (2.5, 9)	0 (0, 0)	0.0002
Respiratory distress 0 (none) to 10 (life-threatening)	5.3 ± 3.25	0.7 ± 1.08	< 0.0001
Quality of life 0 (excellent) to 10 (worst imaginable)	6.4 ± 1.93	1.6 ± 1.47	< 0.0001
Exercise tolerance 0 (normal) to 10 (sedentary)	5.7 ± 2.94	2.1 ± 1.92	0.0001

* Data presented as mean ± SD or median (inter-quartile range) as appropriate

For Veterinary Use Only. Not For Human Use.
Phone: (650) 327-5000
Fax: (650) 332-1600
E-mail: sales@infinitemedical.com
www.infinitemedical.com

