

**Avoid tubing  
misconnections  
with ENFit® and  
McKesson Brands**

# Enteral connections are changing

## It's time to move to ENFit®

Tubing misconnections threaten the safety of patients with some outcomes leading to death. Enteral misconnections are possible when inadvertent connections occur between enteral feeding systems and non-enteral systems such as IV lines, tracheostomy tubes, peritoneal dialysis catheters etc. The Global Enteral Device Supplier Association's (GEDSA) StayConnected® program is working to help combat these occurrences and make enteral feeding safer with the transition to the ENFit® connector.

ENFit connectors are specifically devised so that they screw together and do not allow connectivity with any other connector unrelated to enteral (tube) feeding. This specialized connector will help prevent misconnections or accidental disconnections.

The enteral products impacted include enteral syringes, extension sets and feeding tubes/bags. Button-style G-tubes will not change except for the extension sets that connect to them. Manufacturers have already begun moving to this new design and we want to help make the transition as easy for you as possible.

### When is this change happening?

Previous designs of feeding tubes and cross-application adaptors stopped being made by most manufacturers July 1, 2021. Most manufacturers are discontinuing producing transition sets and adapters effective January 1, 2022.


### Who will be impacted?

- Skilled nursing facilities
- Home health agencies
- Home medical equipment providers
- Home infusion providers
- All other providers providing enteral nutrition therapy

### What's changing for you?


The enteral feeding products you buy today will be changing. As manufacturers move to the ENFit design, connectors will need to be ENFit-compatible.

#### Old connector



**Feeding Tube**  
(Current)

#### "Step" connector



**Transition Set** (Temporary)  
Allows fitment to current feeding port until new ENFit enteral feeding tubes are available.

#### Final connector



**Feeding Tube** (Final)  
Changing from male – the stepped or Christmas tree connector – to the new ENFit female connector. The feeding tube port for the administration set will change from female to male.

Source: Used with permission from Global Enteral Device Supplier Association's (GEDSA)

**Manufacturers and McKesson Brands, have moved to this new design and we want to help make the transition as easy for you as possible.**



## McKesson Brand Assortment of ENFit® products

McKesson Brands offers an assortment of enteral feeding supplies featuring ENFit® connectors. These specially designed tips meet ISO 80369-3 compliancy design standards, supporting patient safety by eliminating risks of IV misconnections between enteral syringes, extension sets and feeding tubes/bags.

Syringes				
Description	MFR #	Item #	QTY	
McKesson Piston Irrigation Syringe with EnFit® Tip, 60 mL, Sterile	911	988343	1/ea, 50/cs	
McKesson Piston Irrigation Syringe, Flat Top, with ENFit® Tip, with IV Pole Bag, Non-Sterile	910	988342	1/ea, 30/cs	

Gravity Feeding Bags				
Description	MFR #	Item #	QTY	
McKesson Enteral Feeding Gravity Set, Small Bore, 1200 mL	16-702500	1109552	1/ea, 30/cs	
McKesson Enteral Feeding Gravity Set, Large Bore, 1200 mL	16-702505	1109551	1/ea, 30/cs	

Extension Sets				
Description	MFR #	Item #	QTY	
McKesson Extension Feeding Set, Dual ENFit® Port with Clamp, 12", Right Angle, Continuous Feeding	194-0121-12	1183126	1/ea, 5/bx, 6 bx/cs	
McKesson Extension Feeding Set, Dual ENFit® Port with Clamp, 24", Right Angle, Continuous Feeding	194-0121-24	1183127	1/ea, 5/bx, 6 bx/cs	
McKesson Extension Feeding Set, Single ENFit® Port with Clamp, 12", Straight, Bolus Feeding	194-0123-12	1183124	1/ea, 5/bx, 6 bx/cs	
McKesson Extension Feeding Set, Single ENFit® Port with Clamp, 24", Straight, Bolus Feeding	194-0123-24	1183125	1/ea, 5/bx, 6 bx/cs	
McKesson Extension Feeding Set, Single ENFit® Port with Clamp, 12", Right Angle, Bolus Feeding	194-0124-12	1183122	1/ea, 5/bx, 6 bx/cs	
McKesson Extension Feeding Set, Single ENFit® Port with Clamp, 24", Right Angle, Bolus Feeding	194-0124-24	1183123	1/ea, 5/bx, 6 bx/cs	



**For more information, please contact your McKesson Account Manager or visit us at [mms.mckesson.com/enfit](https://mms.mckesson.com/enfit).**

## Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

IMPORTANT NOTE: All material provided is intended for informational purposes only and should not be used to replace regulatory or company-specific documents, nor should it replace the advice of a qualified professional. This information is based on the proposed ISO 80369 standard and does not imply or suggest any regulatory clearance of any specific product and may not be predictive of clinical outcomes.