

VLU Care Made Simple



Control for Every Venous Leg Ulcer

VLU Treatment

1. Cleanse

Never compromise **efficacy** or **safety** and take control of the wound environment.

vashe[®] wound solution

Vashe[®] Wound Solution, powered by a proprietary electrochemical process called **Natural Zone Technology 5.5**, is titrated to a skin-friendly pH maximum of 5.5 to mimic the body's natural antimicrobial defenses and provide an ideal healing environment. Formulated with your patient in mind, Vashe remains non-cytotoxic unlike solutions, such as Dakin's solution, that can contain sodium hypochlorite which is known to be cytotoxic even at small concentrations.^{1.4}

Set the ideal stage for endogenous healing by using Vashe first.

Never compromise.



2. Dress

Fight **moisture** and **infection** in two ways.



Drawtex Hydroconductive Wound Dressing used with either UrgoTul or UrgoTul Ag can accomplish:

1. Active removal of excess moisture and impediments that compromise wound healing from the deep tissue.⁵⁻¹⁰

2. Flexibility to offer a moisture balanced antibacterial barrier on the wound surface to reduce infection.¹¹

Take control of the wound environment with an active, adaptable dressing solution.



3. Compress

Compression application is not the time for guesswork.

UrgoK2[™] & UrgoK2[™]

UrgoK2 Dual Compression System removes the variables of bandage tension, bandage overlap, and limb circumference to ensure compression that is:

Consistent: 87% of nurses achieved therapeutic compression on their very first application¹²

Continuous: Short stretch and long stretch bandages uniquely combine to provide therapeutic compression day and night, whether the patient is in motion or not¹³

 $\label{eq:comfortable: 95\% of patients preferred UrgoK2 during the day and 92\% during the night^{14}$

Trust that your compression system is spot on.





Urgo's select portfolio of industry-changing products puts you in control of your patients' VLUs.

Urgo Product Information

UrgoK2[™] Not made with natural latex

Size	Width	Dosage	Code
7⅓-9¾ in (18-25 cm)	4 in (10 cm)	40mm/Hg	553243
9¾-12½ in (25-32 cm)	4 in (10 cm)	40mm/Hg	553244

Size	Width	Dosage	Code
71⁄8-93⁄4 in (18-25 cm)	4 in (10 cm)	20mm/Hg	553245
9¾-12½ in (25-32 cm)	4 in (10 cm)	20mm/Hg	553246

UrgoTul™

Size	No. per box	Code
2 x 2 in (5 x 5 cm)	10	506487
4 x 5 in (10 x 12 cm)	10	506488
6 x 8 in (15 x 20 cm)	10	506489

UrgoTul[™]Ag/Silver

	Size	No. per box	Code	
	4 x 5 in (10 x 12 cm)	10	509341	
	6 x 8 in (15 x 20 cm)	10	509342	
	8 x 16 in (20 x 40 cm)	5	553228	

drawtex A Hydroconductive Wound



* Tracheostomy and Tube Dressing

vashe

wound solution			
Bottle Size/Pack Size	Quantity	Code	
4.0 fl. oz. (118 mL)	Bottles/24-Pack	00312	
8.5 fl. oz. (250 mL)	Bottles/12-Pack	00313	
16.0 fl. oz. (475 mL)	Bottles/12-Pack	00314	
34.0 fl. oz. (1 liter)	Bottles/6-Pack	00322	

References: 1. Wang L, Bassiri M, Najafi R, et al: Hypochlorous acid as a potential wound care agent: Part I. Stabilized hypochlorous acid: A scomponent of the inorganic armamentarium of innate immunity. J Burns Wounds. 2007;6:65-79. 2. Block SS. Disinfection, sterilization, and preservation. Philadelphia: Lea & Febiger; 2000. 3. Sampson CM, Sampson MN. Hypochlorous acid: A safe and efficacious new wound therapy. Poster presented at: World Union of Wound Healing Societies; 2008; Toronto, Ontario, Canada. 4. Hidalgo E, Bartolome R, Dominguez C. Cytotoxicity mechanisms of sodium hypochlorite in cultured human dermal fibroblasts and its bactericidal effectiveness. Chem Biol Interact. 2002;139(3):265-282. 5. Smith DJ, Karlnoski RA, et al. The treatment of partial-thickness burns with a hydroconductive wound dressing: clinical and mechanistic effects. Surgical Science. 2013;4:268-272. 6. Wolvos T, Livingston, M. Wound Fluid Management in Wound Care: The Role of a Hydroconductive Dressing. Wounds. 2013; 25(1): 1-8. **7.** Wolvos T. Analysis of wound bed documentation in advanced wound care using Drawtex, a hydroconductive dressing with LevaFliber technology. Wounds. 2012;24(9)(Suppl):9-10. **8.** Wolcott RD, Cox S. The effects of a hydroconductive dressing on wound biofilm. Wounds. 2012;24(9) (Suppl):14-16. **9.** Ortiz RT, Moffatt LT, et al. In vivo and in vitro evaluation of the properties of Drawtex LevaFliber wound dressing in an infected burn wound model. Wounds. 2012;24(9) (Suppl):14-13. **11.** Lazareth I., Meaume S., Sigal-Grinberg M-L, et al. The Role of a Silver Releasing Lipido-colloid Contact Layer in Venous Leg Ulcers Presenting Inflammatory Signs Suggesting Heavy Bacterial Colonization: Results of a Randomized Controlled Study. WOUNDS 2008;20(6):158–166. **12.** Pilati, L. Comparing 2-Layer with Traditional 4-Layer Compression Therapy. Ostomy Wounds Management. 2002; 66(4): **8** & 10. **13.** Jünger M, et al. Comparison of interface pressures of three compression bandaging system sues on healthy volunteers. J Wound

Distributed by: Urgo Medical North America, LLC. 3801 Hulen Street, Suite 251, Fort Worth, TX 76107 To order, call 1-855-888-8273 or visit www.urgomedical.us

© 2021 Urgo Medical North America. All rights reserved. Vashe, UrgoTul, UrgoTul Ag, UrgoK2, Urgo, and the Urgo logo are registered trademarks of Urgo Medical. Drawtex is a registered trademark of Beier Drawtex Healthcare (Ptv) Ltd.

COLLTOO4 Rev. 08/21