

An innovative benchtop platform for automation and miniaturization of a PCR-free WGS library preparation protocol for Salmonella enterica NGS typing

Sebastian Aguilar Pierlé¹, Ouerdia Maskri¹, Camille Soucies¹, David Gosselin¹, Amel Bendali¹ ¹Inorevia, Paris 75014

KEY BENEFITS

INTRODUCTION

Salmonella enterica is the second most reported bacterial cause of food-borne infections in Europe. Therefore, surveillance activities based on pathogen subtyping are important for controlling Salmonellosis. Whole-genome-based typing is becoming the gold standard due to its improved resolution and reduced cost.

Here, we showcase the advantages of using Magelia's PCR-free library preparation application for blind NGS typing of a relevant foodborne pathogen. Our multi-OMICs platform challenges the minimum required input for a PCR-free library kit, by reducing reagent volumes 8x and DNA input 50x while yielding high quality sequencing libraries

PCR-FREE WGS LIBRARY PREP USING STANDARD HUMAN gDNA

PCR-free WGS library preparation was performed using the Kapa Hyper Prep PCR-free kit (Roche) (Figure 1). Standard human gDNA from Promega was treated at three different input concentrations following the manufacturer's instructions and using in parallel three Magelia prototypes (Table 1)



Salmonella enterica NGS TYPING APPLICATION

In order to validate our PCR-free application with clinical field samples, eight different S. enterica serotypes were used (Table 2). We chose the two most widely prevalent serotypes (Typhimurium and Enteritidis), as well as closely related ones to highlight the resolution of the genomic data generated.

Library preparation was performed as previously described and samples were treated blindly, in parallel in the Magelia and manually (Figure 1). Due to the clinical origin of the samples, an assortment of inputs was treated (Table 2).

Prepared libraries underwent standard Quality Control. Samples were then normalized, pooled then sequenced on a MiSeq[™] Sequencing System with a 2x75 bp configuration.

Salmonella NGS TYPING PRE-SEQUENCING METRICS

Magelia treatment of S. enterica samples resulted in higher library yield compared to manual preparation for all but one sample, sample U (Figure 3A). This was due to the fact that PCR-free treatment for this manually processed sample did not yield enough material for sequencing, hence, amplification was necessary.

Remarkably, average adapter dimer percentages were significantly lower (p=0,04) in Magelia treated samples when compared to manual treatment (Figure 3B).

S. enterica NGS TYPING AND POST-SEQUENCING METRICS

We showed that 90,39 % of reads for the manually treated samples and 95,67 % of the reads for the Magelia treated samples have a PHRED score above 30, implying a larger proportion of high-quality reads (Figure 4). Moreover, the manually treated samples display a more spread distribution for reads with scores between 20 and 38, confirming lower overall quality for these samples.

Higher quality sequencing data resulted in improved coverage and assemblies for Magelia treated samples. This in turn, allowed for detection of more genetic elements of clinical relevance (Figure 5).

An epidemiologically relevant SNP was detected for the Panama isolate in the Magelia treated sample and not in its manual counterpart. This was due to a lack of coverage in the region of interest for the manually treated sample (Figure 6).

Suitable coverage obtained for Magelia treated samples allowed us to successfully identify all serotypes. Conversely, only 5/8 were identified through manual treatment (Table 3).

Identifier	Magelia	Manual
U	+	+
S	+	+
Y	+	+/-
w	+	-
н	+	+
N	+	+
Р	+	+
R	+	-
	U S Y W H N P	U + S + Y + M + H + N + P +

underline a higher success rate (relevant for regulated environments)

CONCLUSIONS

nucleotide resolution:

technology

High Quality data from Reagent volume low input samples reduction Full automation of complex workflows ğ Compact benchtop platform Turn-around time and High Efficiency hands-on-time reduction reactions Automated reagent dispensing 1 困 Ì

> Prepared libraries were quality checked using the Bioanalyzer® and Qubit®. As shown in Table 1, Magelia treated samples yielded 3x more ready-to-sequence libraries when compared to manual treatment.


	Manual		Magelia		Manual		Magelia		Manual		Magelia							
Sheared DNA Input in ng	50ng	50ng	50ng	50ng	50ng	50ng	20ng	20ng	20ng	20ng	20ng	20ng	5 ng	5 ng	5 ng	5 ng	5 ng	5 ng
Magelia prototype				Magelia 1	Magelia 2	Magelia 1				Magelia 2	Magelia 1	Magelia 3				Magelia 2	Magelia 3	Magelia 1
Final Concentration in nM	15,8 nM	12,9 nM	12,2 nM	51,5 nM	52,4 nM	51,9 nM	6,1 nM	6,4 nM	5,6 nM	18,9 nM	25,4 nM	19,7 nM	1,45 nM	1,3 nM	n/a	3,16 nM	5,1 nM	4,4 nM

Table 1. WGS PCR-free Library preparation using standard human gDNA

Importantly, using as little as 5 ng as input (10x less than the kit's recommended limit), yielded sufficient material for sequencing.

Furthermore, final library concentrations were comparable between the three prototypes, highlighting Magelia's robustness.

Serovar	Identifier	DSMZ nb	Isolation	Origin	Input qty				
Typhimurium	U	101475	Copper fed pigs 2013	Denmark	5 ng				
Gallinarum	S	4883	?	?	35 ng				
Enteritidis	Y	107306	Patient with diarrhea	Georgia	35 ng				
Saintpaul 4,5:e,h:1,2	w	27656	?	?	35 ng				
Agona	н	102864	?	?	85 ng				
Senftenberg (N)	Ν	10062	?	?	35 ng				
Heidelberg	Р	9379	?	Germany	35 ng				
Panama (1,9,12:1,v:1,5)	R	9145	Baby 1994	USA	35 ng				
Table 2. Salmonella serovars treated blindly for NGS typing									



genomic region containing SNP c.884A>C

MAGELIA, A UNIQUE multi-OMICs PLATFORM:

Magelia treatment of low input samples enabled characterization of clinically relevant pathogens to the single

- Disruptive core technology enabling high precision magnetic bead handling No evaporation, allowing for improved reaction
 - kinetics
- Optimal cleanup and negligible adapter-dimer formation made possible thanks to our patented magnetic tweezer Higher amplification efficiency: less PCR cycles needed, reducing bias
- Perspectives: we are continuously developing a variety of miniaturized applications for low input samples: o RNA-seq library prep with rRNA depletion/mRNA capture o Chromatin Immuno-Precipitation (ChIP) followed by NGS library prep

• Blind identification of all Salmonella serotypes for Magelia treated samples vs. 5/8 samples for manual treatment

Faster reactions (2-hour immunoprecipitation in Magelia vs. 16 hours for manual preparation)

PLATFORM FEATURES

- Integrated heating/cooling/thermal cycling

- User defined protocol parameters
- Proprietary cartridges prevent inter and

intra run cross-contamination