



TPASS Visitor Management System – Enhancements to support your efforts against the Coronavirus

.....

Thermal Imaging Camera Platforms

**Eyemetric Identity Systems
2020**



New Health Check Status Features

The screenshot displays the Eyemetric Identity Systems Visitor Management interface. The top navigation bar includes links for Home, Student, Visitor, Personnel, Reports, and Administration. The 'Visitor' tab is active. A search bar contains the text 'boll'. To the right of the search bar are buttons for '+ Host', 'Edit Visitor', and '+ New Visitor'. The main form area shows a visitor's profile for 'RAYMOND L BOLLING III'. The visitor's birth date is '6/10/1965' and their type is 'VISITOR'. A red box highlights the visitor's photo, which is labeled 'Flagged' with a '2' in a yellow box. A yellow box with the number '1' highlights the 'Health Check' section, which has radio buttons for 'Passed' and 'Failed' (selected), a date of '8/4/2020', and a remark 'Family Member Tested Positive'. A yellow box with the number '3' highlights a green box at the bottom left that says '0 of 40 Max'. At the bottom of the form are two buttons: 'Sign-In Visitor' and 'Cancel'.

eyemetric[®]
IDENTITY SYSTEMS

Home Student Visitor Personnel Reports Administration

Visitor Student Visitor

Search boll

+ Host Edit Visitor + New Visitor

First name Middle name Last name

RAYMOND L BOLLING III

Birth date 6/10/1965 Type VISITOR

Remarks

Location Main Office

Purpose Select...

Health Check ☐ Passed ☒ Failed

Date 8/4/2020

Remarks Family Member Tested Positive

Flagged

0 of 40 Max



Sign-In Visitor Cancel

1. Added the ability to record visitor's health screening status.
2. Able to "Flag" a visitor's barring them access for a set number of days
3. Added the ability to create a maximum number of visitor's allowed in at one time



Integration with Thermal Body Cameras

Integrating Thermal Screening “Seek Scan” with TPASS VMS and Face Recognition



Visitor

Student Visitor

Q Name or ID of the Visitor...

Host Edit Visitor New Visitor

First name Middle name Last name

RAYMOND L BOLLING III

Birth date 6/10/1965 Type VISITOR


Remarks

Location Select...

Purpose Select...

Sign-In Visitor Cancel

Photo	Name	Relationship	Address
-------	------	--------------	---------



Integration with Thermal Body Cameras



Fotric 226B Camera



TPASS VMS processes the thumbnail picture against FR database for a match and automated alert being sent out by TPASS VMS.



Integration with Thermal Body Cameras

Fotric 226B Camera System – Mobile Screening Cart




Fotric offers a complete bundle with large display, laptop, rolling stand and secondary thermometer.



Integration with Thermal Body Cameras

Integrating Thermal Screening Kiosk with TPASS VMS and Face Recognition



eyemetric
IDENTITY SYSTEMS

Home Student Visitor Personnel Reports Administration

Visitor Student **Visitor** Q doe X + Host Edit Visitor New Visitor

First name Middle name Last name
JANE DOE

Birth date Type PARENT

Remarks

Health Check ☒ Yes ☐ No

Date 7/7/2020

Remarks Family Member has tested positive for Coronavirus

Location Main Office

Purpose Select...

0 of 50 Max

Sign-In Visitor Cancel

