



Improved Mobility. **Made Easier.™**

L300 Go is a functional electrical stimulation (FES) system that satisfies the productivity demands of today's value-based healthcare system. Key aspects of the L300 experience have been dramatically improved with 3D Motion Detection, multi-channel stimulation, Bluetooth® programming and a home user mobile app that tracks activity to keep patients engaged in the rehabilitation process. All of this in a streamlined design, with a fitting process that is faster and easier than ever before.



3D MOTION DETECTION



MULTI-CHANNEL STIMULATION



FAST, INTUITIVE SET-UP



OUTCOME MEASURES



HOME USER APP

IMPROVED CLINIC EXPERIENCE!

NEW CUFF DESIGN

15% WEIGHT REDUCTION

IMPROVED COMFORT



LAURI L.
Stroke Survivor



It's **Go** Time!

Technology has advanced and the demands being placed on healthcare providers have changed. Five revolutionary changes make **L300 Go** more capable and versatile.



3D MOTION DETECTION

L300 Go is the first FES system to monitor movement in all three planes using a 3-axis gyroscope and accelerometer. A new learning algorithm adapts to changes in gait dynamics and deploys stimulation in 0.01 seconds. **That's reliable response - 30 times faster than the blink of an eye.**



MULTI-CHANNEL STIMULATION

Precisely control the amount of dorsiflexion and inversion/eversion - with a single electrode. **It's easier than ever to optimize gait dynamics - no 'fishing' required!**



FAST, INTUITIVE SET-UP

Integrated motion detection and onboard controls make the Foot Sensor and Control Unit optional. Fewer components to manage, a Quick Start Fitting Mode and Bluetooth® programming **significantly reduce set-up time and maximize the productivity of each therapy session.**



OUTCOME MEASURES

Use the integrated 10 Meter Walk Test to easily track and document the progress needed to **succeed in today's value-based healthcare model.**



HOME USER APP

The **myBioness™** mobile app allows users to set goals, monitor activity and see how far they've come, **allowing the gains made in therapy to continue.**

CONTACT US AT 800.211.9136, OPTION 2 | [L300Go.com](https://www.l300go.com)

Individual results vary. Consult with a qualified physician to determine if this product is right for you.

Important Safety Information and Risks: For Indications for Use, Contraindications, Warnings, Adverse Reactions, Precautions, and other safety information please refer to www.bioness.com/Safety_and_Risk_Information.php (also available in the L300 Go Clinician's Guide).

L300 Go®, myBioness™, Bioness, the Bioness Logo®, LiveOn® and Improved Mobility. Made Easier.™ are trademarks of Bioness Inc. in the United States or other countries | www.bioness.com | **Rx Only**



Bioness Inc.
25103 Rye Canyon Loop
Valencia, CA 91355 USA
Telephone: 800.211.9136 or 661.362.4850
Website: www.bioness.com



Bioness Europe B.V.
Stationsweg 41
3331 LR Zwijndrecht, The Netherlands
Telephone: +31.78.625.6088
Email: international@nl.bioness.com
Website: www.bioness.com