

Since our inception in 1998, A. Titan has been recognized as an innovative leader in the dental community. It's our mission to enhance the quality of patient care by providing superior quality products and instrumentation on the cutting edge of dentistry.

Focused on expanding the reach of the A. Titan brand, our diverse product line is offered through dealers around the world. Headquartered in Orchard Park, New York, our instruments are produced in our 100,000 square foot state of the art manufacturing facility or designated for production in Germany at our partner location.

We strive to achieve our mission of enhancing patient care by providing continuing education support to universities throughout the world, working closely with key opinion leaders to develop innovative instrumentation, and maintaining the highest levels of customer service. Forging long lasting relationships, our team focuses on bringing products to market that encourage the promotion of aesthetic perfection.



Banditt X-TRA GRIP Forceps

Featuring a Diamond NANO Composite Coating



150XG Upper Universal Forceps



151XG Lower Universal Forceps



7990XG Lower Universal Forceps



7301XG Lower Molar Ash Pattern Forceps



1XG Upper Anterior Forceps



3608XG Lower Anterior Forceps

- Unique properties of natural diamonds
- High micro hardness
- Exceptional corrosion resistance
- Abrasive wear resistance for up to 50 years
- Reduces hard and soft tissue damage
- Beak innovation with biomechanical grooves
- Ergonomic, lightweight design reduces hand fatigue
- Recessed lateral serrations provided an enhanced grip
- Made in Germany

PLUS SERIES FORCEPS





51+ Upper Root Forceps



1+ Upper Anterior Forceps



150+ Upper Universal Forceps



151+ Lower Universal Forceps



210+ Upper Molar Forceps

- Lightweight and balanced with an adult handle ideal for primary teeth
- Beveled and serrated beaks allow for more tooth and root surface contact and operator control
- The most advanced design using CAD CAM technology and experienced Board Certified Consultants
- Made from 100% stainless steel
- Crafted in Germany

RAC Forceps

EXTRACTION ATRAUMATIC

Exclusive Handle Designs:

- Make forceps 25% lighter
- Offer great ergonomic comfort
- · Provide increased tactile sensitivity

Unique Beak Designs:

- Preserve bone
- · Provide an exceptional grip of the tooth



Upper Universal Forceps 0700



Upper Anterior Forceps-Notched 3400



Upper Universal Forceps-Notched 3500



Lower Anterior Forceps-Notched



Lower Universal Forceps-Notched 2190



3608



Upper Universal Forceps 3508



Lower Universal Forceps



Upper Molar Left Forceps 1700



Lower Cowhorn Forceps





Upper 3rd Molar Forceps

Upper Molar Right Forceps



3-Prong Upper Left Molar Forceps 3-Prong Upper Right Molar Forceps



Lower Molar Forceps 7301

Lower 3rd Molar Forceps



Universal Maxillary Root Forceps 4415



Upper Root Forceps 5115



Lower Root Forceps 4515





Lower Narrow Root Fragment 4507



Cowhorn, Lower Molar Forceps



ATRAUMATIC EXTRACTIONS

Beveled tips follow root structure down to PDL. Using 2-3 degree rotation around circumference of tooth, you will gently separate the fibers from the root structure down to the PDL, creating a trough that will minimize tissue and bone trauma.

XO-10

Small straight tip works on anterior teeth

XO-20

Small curved tip works on anterior and single rooted teeth

XO-30

Small angled tips work on the mesial of anteriors & bicuspids

XO-40

Small angled tips work on the distal of anteriors & bicuspids

XO-50

Medium angled tips work on 1st, 2nd, and Premolars

XO-60

Medium angled tips work on 1st, 2nd, and Premolars

XO-80 Straight Micro

XO-110

- Hybrid instrumentation, a cross between an Elevator and a Periotome
- Custom sure-grip handle for controlled luxation force
- Honed narrow tips to disrupt the periodontal fibers
- Preserves tissue structures
- Eliminates bone fracture and promotes minimally invasive exodontia
- · Made from 100% stainless steel
- Made in America

XO-70

Straight spade shaped blade ideal for interproximal tissue reflecting

XO-77

Beveled spade shaped blade slides nicely below gum line and follows the contour of root to both sever fibers and aid in root removal

XO-77D

Angled for the distal, the beveled spade shaped blade slides nicely below gum line and follows the contour of root to both sever fibers and chase ankylosed roots

XO-77M

blade slides nicely below gum line gled for the mesial, the beveled spade shaped blade slides nicely below gum I and follows the contour of root to both sever fibers and chase ankylosed roots

XO-77AL

Angled left, the beveled spade shaped blade slides below gum line to sever fibers and chase ankylosed roots

XO-77AR

Angled right, the beveled spade shaped blade slides below gum line to sever fibers and chase ankylosed roots


XO-90

Large angled tip excellent for 3rd molars

XO-100

Large angled tip gives the ability to maneuver and elevate tuberosity around 3rd molar





Elemental Kit K-3300

Contains:

XO-10, XO-20, XO-77 and a sterilization cassette

X-Otome Kit K-6600

Includes one of each: XO-10, XO-20, XO-30, XO-40, XO-50, XO-60 XO-77 and a sterilization cassette

X-Otome Kit K-7700

Includes one of each: XO-77, XO-77AL, XO-77AR, XO-77M, XO-77D and a sterilization cassette

BONE PRESERVATION ELEVATORS

- Minimize bone loss
- Create space between the bone and the root for all atraumatic extractions
- Thin sharp tip of the elevator slides between the bone and the root creating space without removing any bone
- Fast and easy removal of root and fragments without a surgical procedure
- Made in Germany from 100% stainless steel
- Lifetime Guarantee







Bone Preservation Elevator Kit K-5000

Contains:

#63, #60B, #60R and #60L

RIGHT ANGLED DRIVER

Our Right Angled Driver makes it easy to place posterior screws in the mandible or maxilla. Because it's manually driven, there's no need to connect it to a hand piece.



- Maximize access and comfort with the ability to rotate the head to the exact angle needed
- · Make quick and easy abutment changes in posterior
- Manually driven, no need to connect to hand piece



ATA Torque Adapter

- 7 torque values 10, 15, 20, 25, 30, 32, 35 Ncm
- Easily connects to Right Angled Driver (RAD)
- Never over torque. RAD makes a "clicking" sound when desired Ncm is achieved
- Easily set Ncm to place desired implant
- Autoclavable



BITS KITBits for Right-Angled Drivers,
Set of 7

- Compatible with all major implant systems
- Multiple bit options cover all existing and new patients
- Bits are used with Right Angled Driver (RAD) or any standard latch-type contra-angle system
- Ideal for posterior abutments and cover screws in the mandible and maxilla
- Each bit is laser etched for quick identification

SCISSORS

- Super Cut Scissors feature one serrated blade which provides a smooth cut and stay sharper longer
- Excellent for trimming tissue and suture material
- Made in Germany from 100% stainless steel







List of Scissors

808SC - Goldman Fox, 5"

840SC - Iris Curved, Extra Fine 4.5"

865SC - Iris Curved, 4.5"

869SC - Iris Straight, 4.5"

877SC - LaGrange, 4.5"

8111C - Iris Curved, Carbide Tip, 4.5"

8110C - Iris Straight, Carbide Tip, 4.5"

840 - Iris Curved, Extra Fine, 4.5"

865 - Iris Curved, 4.5"

869 - Iris Straight, 4.5"

808 - Goldman Fox, Serrated, 5"

818 - Metzenbaum 5.75"

815 - Rounded Tip, 4.5"

877 - Curved LaGrange, 4.5"

859 - Curved Kelly, Serrated, 6.25"

861 - Dean Scissor, Serrated, 7"

880 - Closed Handle Locklin, Serrated, 6.25"

882 - Open Handle Locklin, Serrated, 6.25"

899 - Straight Suture Scissor w/ Loop, 6"


809 - Curved Suture Scissor w/ Loop, 4.5"

378 - Crown & Bridge, 4"

878 - Crown & Bridge, 4.5"

8546SC - General Office


NEEDLE HOLDERS



- Fully Autoclavable
- 100% stainless steel
- Made in Germany

6"





5.5"

List of Needle Holders

M7 - Mayohegar, 7"

M6 - Mayohegar, 6"

Micro 6C - Micro, 6"

C6C - Crilewood, 6"

B5 - Baumgartner, 5.5"

D4 - Derf, 4 3/4"

S6 - Swedish, 6"

06C - Olsen-Hegar, 6.25" (Scissor/Needle Holder Combo)

05C - Olsen-Hegar, 5.5" (Scissor/Needle Holder Combo)

MA14C - Mathieu, 5.5"

MR17C - Mathieu-Ryder, 6.75"

340TC - Boynton, 4.75"

188TC - Flat Handle, 5.5" Castroviejo

188TC-LONG - Flat Handle, 7" Castroviejo

188TC-CURVED - Flat Handle, 5.5" Castroviejo - Curved Beaks

2000 TC - Titan Designed Handle, 5.9" Castroviejo

187TC - Round Handle, 5.5" Castroviejo

187TC-LONG - Round Handle, 7" Castroviejo

187TC-CURVED - Round Handle, 5.5" Castroviejo - Curved Beaks

RS3 RAZOR SHARP SUTURES

- Smooth, polished needle, easy glide surface, and uniform thickness
- Superior suture material for advanced wound closure
- Minimal needle to material transition reduces suture line bleed
- · No coagulation of blood on suture
- 12 sterile sutures per box
- Made in America




						Total Co.
	A. TITAN SUTURE CODE	NEEDLE	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE	ETHICON CODE
ABSORBABLE PGA Violet & Braided	ATI-6-18-3PGAV	C-6	3/8 Circle 19mm	18"	3-0	J393H
	ATI-6-18-4PGAV	C-6	3/8 Circle 19mm	18"	4-0	J392H
	ATI-3-18-5PGAV	C-3	3/8 Circle 13mm	18"	5-0	J463G
	ATI-6-18-5PGAV	C-6	3/8 Circle 19mm	18"	5-0	J391H
ABSORBABLE						
PGA Fast Absorbing	ATI-31-18-3PGAR	C-31	3/8 Circle 19mm	18"	3-0	VR497
& Undyed	ATI-31-18-4PGAR	C-31	3/8 Circle 19mm	18"	4-0	VR496
a onayou						
	ATI-6-27-3PGA	C-6	3/8 Circle 19mm	27"	3-0	J423H
	ATI-31-18-3PGA	C-31	3/8 Circle 19mm	18"	3-0	J497G
	ATI-6-27-4PGA	C-6	3/8 Circle 19mm	27"	4-0	J422H
ABSORBABLE	ATI-3-18-4PGA	C-3	3/8 Circle 13mm	18"	4-0	J494G
PGA Undyed	ATI-31-18-4PGA	C-31	3/8 Circle 19mm	18"	4-0	J496G
& Braided	ATI-3-18-5PGA	C-3	3/8 Circle 13mm	18"	5-0	J493G
	ATI-31-18-5PGA	C-31	3/8 Circle 19mm	18"	5-0	J495G
	ATI-3-18-6PGA	C-3	3/8 Circle 13mm	18"	6-0	J492G
	ATI-3-18-4CG	C-3	3/8 Circle 13mm	18"	4-0	1654G
ABSORBABLE	ATI-3-18-5CG	C-3	3/8 Circle 13mm	18"	5-0	687G
Chromic Gut	ATI-6-27-3CG	C-6	3/8 Circle 19mm	27"	3-0	636H
Officialic dut	ATI-6-27-4CG	C-6	3/8 Circle 19mm	27"	4-0	635H
	711021100		0,0 0,00 10,00			00011
	ATI-6-27-3PG	C-6	3/8 Circle 19mm	27"	3-0	H822H
	ATI-6-18-4PG	C-6	3/8 Circle 19mm	18"	4-0	N/A
ABSORBABLE	ATI-3-18-5PG	C-3	3/8 Circle 13mm	18"	5-0	686G
Plain Gut	ATI-6-27-4PG	C-6	3/8 Circle 19mm	27"	4-0	H8218
	ATI-7-27-3PG	C-7	3/8 Circle 24mm	27"	3-0	H810H
	ATI-9-18-3PG	C-9	1/2 Circle 23mm	18"	3-0	612G
NON-ABSORBABLE Polypropylene	ATI-1-18-6PPN	C-1	3/8 Circle 11mm	18"	6-0	8697G
	ATI-1-18-7PPN	C-1	3/8 Circle 11mm	18"	7-0	8696G
	ATI-3-18-5PPN	C-3	3/8 Circle 13mm	18"	5-0	8698G
	ATI-3-18-6PPN	C-3	3/8 Circle 13mm	18"	6-0	8695G
	ATL C 40 0DDC		0/0 Cinal - 10	10"	0.0	NI/A
NON-ABSORBABLE Black Braided Silk	ATI-6-18-3BBS ATI-6-18-4BBS	C-6 C-6	3/8 Circle 19mm 3/8 Circle 19mm	18" 18"	3-0 4-0	N/A 683G
	ATI-9-18-3BBS	C-6 C-9	1/2 Circle 19mm	18"	3-0	632G
			1/2 Officie 23/11/11			0320
	Precision Revers	e Cutting		Reverse	Cutting	

11mm 3/8 Circle 13mm 3/8 Circle C-31 PS-2 19mm 3/8 Circle

19mm 3/8 Circle 23mm 1/2 Circle 24mm 3/8 Circle

COMPOSITE INSTRUMENTS



SYRINGES

Aspirating syringes have ergonomic thumb rings that are available in small, medium, and large sizes to enhance comfort and maximize control. The lightweight design reduces hand stress and has no removable parts.

- Comfort grip finger rest and thumb ring
- Made in America

HOW IT WORKS

The threaded tip, barrel, rod/plunger has a harpoon attached that embeds itself into the rubber stopper of the carpule. The harpoon allows the dentist to aspirate (draw back) the injection site to ensure the needle tip is in the correct location before injecting the anesthetic solution.



ASPR-1 Aspirating Syringe

- Holds 1.8cc
- Lightweight comfort fitting ring that provides a secure grip for multiple hand sizes

HOW IT WORKS

The self-aspirating mechanism works by a little bump at the end of the syringe compressing the diaphragm of the anesthetic cartridge during injection. When the pressure of injection is released there is the relaxation of the diaphragm, thus creating negative pressure.

With a self-aspirating syringe, aspiration can occur only after the injection has been started.



Syringe

LIFE STEEL Scalers & Curettes

Tips feature **Ion Charged Technology** to increase edge retention

- #6 Handle
- Available single and double ended
- 440A stainless steel
- · Diamond knurl for optimal grip
- · Able to be re-tipped

- Innovative tip designs
- Guaranteed sharpness
- Computerized hydraulic forming ensures 100% angle accuracy
- · Made in America

GRACEYS*

GR 1-2	GR 11-14
GR 3-4	GR 12-13
GR 5-6	GR 13-14
GR 7-8	GR 15-16
GR 9-10	GR 17-18
GR 11-12	

*All Gracey Curettes available in mini and after 5 patterns

AREA SPECIFIC CURETTES

SUB 0 SUB 13-14 SUB 11-12

NEVI SCALERS

NEVI 1 NEVI 3 NEVI 2 NEVI 4

YOUNGER-GOOD CURETTE

YG 7-8

LANGER CURETTES

Langer 1-2 Langer 5-6 Langer 3-4 Langer 17-18

McCALL CURETTES

McCall 13-14 McCall 17-18 McCall 13S-14S McCall 17S-18S

COLUMBIA CURETTES

COL 2R-2L COL 13-14 COL 4R-4L

BARNHART CURETTES

BARN 1-2 BARN 5-6

SCALERS

N 135 204S N 137 204SD 204 MTJ

SICKLE SCALERS

H5 U 15 U 15-30 H5-U15 H5-J1S H5-47 H5-33 Taylor 2-3 H6-H7

JACQUETTE SCALERS

JAC 1-2 JAC 2-3 JAC 1S-2S JAC 2S-3S

GOLDMAN FOX SCALERS

GF - 5 GF - 6

SUPPORT

ORDERING

- Customer Service Support available Mon-Fri 8:30am 5:00pm EST 1-877-ATITAN1 | Local/Intl : 1-716-667-9211
- 24 hour fax available 716-667-9212
- Email orders@atitan.com or info@atitan.com
- · School/Institution pricing available.
- All online orders are forwarded onto the desired dealer per the request of the customer.

REDEMPTIONS

- To receive free goods in association with a current promotion, please fax or email a copy of the dealer invoice to 716-667-9212 or orders@atitan.com
- Include customer's name and address if not indicated on invoice, and the list of free goods for redemption.
- Promotional goods take 6-8 weeks for processing and delivery.

REPAIRS

- Sharpening and repair service is available for all items, please call customer service for pricing.
- 2-3 week turnaround.
- · Customer incurs shipping cost.
- To start a repair, visit the support section of our website to fill out a form and download a shipping label.

RETURNS

- Instrumentation purchased from A. Titan Instruments can be returned for exchange/credit towards a future purchase order.
- All returns from customers must be in original packaging, unused, and have not been put through the ultrasonic or autoclave. The instruments will be inspected upon return and based on inspection, it will be A. Titan's option to allow the return credited in full or denied. Customer incurs shipping cost(s). All returns are subject to a 15% restocking fee based on received condition of said instrumentation.
- No cash refunds. No returns on custom orders.
- To start a return, visit the support section of our website to fill out a form and download a return label.

WARRANTIES

- Approved products that are deemed defective will be repaired or replaced.
- Warranty doesn't apply to normal wear and tear and for use of instruments beyond their purpose or limitations.
- All Tungsten Carbide needle holders are guaranteed for the first 3 years of normal usage.
- Any modifications or alterations made to instruments other than by A. Titan will render the warranty null and void.
- All warranty claims must be made directly with A. Titan, not the distributor.
- · Warranty doesn't include shipping and handling. Customer incurs all shipping cost.
- Periotome and scaler tips aren't covered under warranty, but can be replaced for a nominal fee.
- Right Angled Driver 1 year warranty internal parts and 3 month warranty on manual latch head.
- Instruments that show wear and tear from normal use are not considered to be defective and are not covered by the warranty.
- Sharpening and repairing minor tip damage are considered routine maintenance and are not covered by warranty.
- Modifying and retipping an instrument or failure to provide proper instrument care, including proper cleaning and maintenance, may void this warranty.





