




BRASSELER USA®

QUARTERLY
PROMOTIONS



JULY-
SEPTEMBER
2021

FINISHING ZIRCONIA CROWNS

Technical Lab Procedure as Taught by Ed McLaren, DDS

SPRUE SEPARATION

859

- For separating sprue from the block, use a thin medium needle diamond
- Available in 5/Pk


8856

- To remove remaining sprue content, use a round-end taper fine diamond
- Available in 5/Pk

Buy 10 Packs, Get 30% Off

Buy 5 Packs, Get 20% Off

Mix and match



Order No.	859	856
Head Size ($\frac{1}{10}$ mm)	018	016
Cutting Length (mm)	10	8.0
31 - FG - Friction Grip (19 mm)		
● Medium	859	—
● Fine	—	8856



SHAPING, TEXTURING & BREAKING GLAZE

Refine occlusal and external anatomy as well as create surface texture with these popular Peter Brasseler Series™ Premium Diamonds and ET® Trimming and Finishing Diamonds.


ET® Trimming & Finishing Diamonds

- Precise needle and bullet shapes great for creating and blending grooves
- Available in 5/Pk


Buy 10 Packs, Get 30% Off

Buy 5 Packs, Get 20% Off

Mix and match



US No.	DET9	DOS2
Order No.	135	247
Head Size ($\frac{1}{10}$ mm)	014	009
Cutting Length (mm)	9.0	3.2
31 - FG - Friction Grip (19 mm)		
● Fine	DET9F	DOS2F
● Extra-Fine	DET9EF	DOS2EF




Peter Brasseler Series™ Diamonds

- Football, Flat-End Taper, Round-End Taper and Flame shapes ideal for smoothing and contouring
- Available in 5/Pk

Buy 10 Packs, Get 30% Off

Buy 5 Packs, Get 20% Off

Mix and match



Order No.	379	847	856	856	862
Head Size ($\frac{1}{10}$ mm)	016	012	012	016	014
Cutting Length (mm)	3.4	8.0	8.0	8.0	8.0
31 - FG - Friction Grip (19 mm)					
● Fine	8379	8847	8856	8856	8862

FINISHING ZIRCONIA CROWNS

Technical Lab Procedure as Taught by Ed McLaren, DDS

FINISHING & POLISHING

Polish all surfaces to a smooth luster.

Dialite® Feather Lite™ Single Layer Polishers

- Ideal for blending of anatomy and to take away sharp edges and scratches left by diamond
- Finger projections help prevent over-polishing

Retail: \$204.50/Pk

**Save 25% When You Buy
3 Or More Packs**

NEW!



Order No.	SFL26CHP*	SFL26MHP	SFL26FHP
Head Size (1/16 mm)	260	260	260
Cutting Length (mm)	1.6	1.6	1.6
Grit	Coarse*	Medium	Fine
Unmounted (Use Mandrel 305.11)			
Qty/Pack	10	10	10

*Use only on post-sintered zirconia

Suggested Tip: Dr. McLaren recommends using at low speed (1,000-1,500 rpm)



Dr. Ed McLaren created the ZR TFT (The Final Touch) Procedure for pre-sintered & sintered zirconia to offer a practical and efficient solution for lab technicians and dental practices with in-office CAD/CAM. The technique, used and taught at the Art Oral America and DTG Education Center, utilizes Brasseler fine diamonds to shape and texture and Feather Lite™ polishers for finishing. It's all the instruments needed to create preparations for indirect restorations for maximum fit and quality. Visit www.EdMcLaren.com to learn more.

NEW! 3111FGHP Converter Mandrel



Retail: \$119.50/Pk


**Save 20% When You Buy
2 Or More Packs**

- ▶ Use with diamonds featured on page 2 to expand your lab instrument options
- ▶ Uses friction grip burs in straight handpiece lab motors (like the Forza L50K Lab Motor)
- ▶ Snap-on clamping sleeve
- ▶ Conveniently bundled 10 per pack so each rotary instrument has its own mandrel
- ▶ Overall length 10mm - significantly shorter than Brasseler USA Converter Mandrel 260BR6082 90 at 44mm
- ▶ Shorter length provides for greater control

3111FGHP Converter Mandrel



260BR6082 90 Converter Mandrel




CROWN PERMANENT CEMENTATION

CLEAN-UP, ADJUSTMENT, & POLISHING

SV3 Savannah-3 Scaler

- Remove excess cement to reduce plaque formation with this specialized instrument.



Save 20% When You Buy 6 Or More Instruments

Mix and match

Retail: \$37.95 ea.

Dialite® Diamond

- Use this specifically designed football-shaped, fine-grit diamond for final adjustments and to minimize fracturing of all modern ceramics.

Buy 10 Packs, Get 30% Off
Buy 5 Packs, Get 20% Off

Mix and match

- Available in 5/Pk



Order No.	369D
Head Size (1/10 mm)	025
Cutting Length (mm)	5.5
31 - FG - Friction Grip (19 mm)	
● Fine	8369DF

Polishing of zirconia to minimize opposing dentition wear can be accomplished by using these popular Brasseler polishers. They're available for all ceramic crown materials including zirconia, lithium disilicate and porcelain.

Dialite® ZR Intra-Oral Feather Lite™ Polishers

Retail: \$62.85/Pk

Special: \$50.28/Pk

Dialite® ZR Intra-Oral Feather Lite™ Assortment Pack

Retail: \$125.70/Pk

Special: \$94.28/Pk

Multi-Use
Two-Step
Polishing System



Order No.	FLW14MZR	FLW14FZR	FLW14AZR
Head Size (1/10 mm)	140	140	140
Thickness (mm)	1.6	1.6	1.6
Grit	Medium	Fine	Medium / Fine
21 - RA - Right-Angle Latch (22 mm)			
Qty/Pack	3	3	Assortment: 6 (3 of each)

Single-Use
Two-Step
Polishing System



Order No.	SFLBCER-1	SFLBF-2	SBFLCP-CE	SBFLCP-F
Head Size (1/10 mm)	100	100	055	055
Grit	Medium	Fine	Medium	Fine
21 - RA - Right-Angle Latch (22 mm)				
Qty/Pack	25	25	25	25

Recommended speed: 4,000-6,000 rpm, Final Polish: 2,000-3,000 rpm

BrioShine™ Intra-Oral Feather Lite™ Polishers

Retail: \$106.95/Pk

Special: \$74.86/Pk

NEW! BrioShine™ Intra-Oral Feather Lite™ Cup Polishers


Retail: \$87.95/Pk

Special: \$61.56/Pk

CROWN REMOVAL

SECTION & REMOVE

Ideal burs and instruments for when the structural integrity of the zirconia cemented crown restoration cannot be preserved or reused.



TKS2 Diamond

- Create a fracture line to section zirconia material quickly and easily with this uniquely designed diamond
- Available in 5/Pk

Buy 10 Packs, Get 30% Off
Buy 5 Packs, Get 20% Off

Mix and match

			
Order No.	TKS2	856	856
Head Size (1/10 mm)	014	016	018
Cutting Length (mm)	8.0	8.0	8.0
31 - FG - Friction Grip (19 mm)			
● Super-Coarse	T35005	—	—
● DuraCut™	—	6856DC	6856DC

DuraCut™ 856 Diamonds

- Unique diamond bonding process
- Highly efficient cutting performance
- Available in 5/Pk


Buy 10 Packs, Get 30% Off
Buy 5 Packs, Get 20% Off

Mix and match

Brasseler has instruments that work to help remove crowns rather than apply a pulling force. Simply place the instrument into separation groove and rotate to spread.

**Save 20% When You Buy
6 Or More Instruments**


Mix and match



NEW! EZ234 Legend Crown Spreader

Retail: \$89.99 ea.

- Ergo Elite Handle
- Multi-angle blades
- Lightweight and ergonomic for better grip



CRZ Zena Crown Spreader

Retail: \$72.25 ea.

- Lightweight and ergonomic
- Non-slip stainless steel handle

INTRODUCING BRASSELER'S HANDPIECE MAINTENANCE SOLUTION...



Designed with convenience in mind, the EverClean Maintenance Machine is the perfect solution for today's busy practice. Combined lubrication and purging cycles lead to a quick and easy, hassle-free maintenance that keeps your handpieces running like new. Affordable, fast, and thorough, the EverClean is the ideal maintenance solution for all your handpiece needs.



Adapters sold separately

- ▶ Complete with three ports for air or electric attachments
- ▶ Dedicated chuck cleaning station to improve bur retention
- ▶ Adapters available for most all handpieces and attachments
- ▶ Minimal clean-up required with internal oil suction feature
- ▶ A fraction of the cost of other maintenance systems

EverClean Maintenance Machine

- 3 cleaning ports
- 1 dedicated chuck cleaning station
- Oil tray and removable door shield
- 2-year warranty

Retail: \$2,199.00

Special: \$1,899.00



EverClean Spray

Retail: \$49.00

Special: \$36.75

Spray can sold separately

ENDOSEQUENCE® FLOWABLE DAM MATERIAL

Introducing EndoSequence Flowable Dam Material, a light-cure, heatless, liquid dam and block-out resin which can be applied to gingiva or tooth surfaces prior to root canal therapy, bleaching, sandblasting, or other procedures requiring intra-oral protection. EndoSequence Flowable Dam Material is provided in pre-loaded syringes and is dispensed using minimum waste precision applicator tips. Upon completion of the procedure, the flexible blue material can quickly and easily be removed with a hand instrument.

- ▶ Flowable yet stackable consistency
- ▶ Light-cures quickly and stays cool to the touch
- ▶ Flexible and tear-resistant
- ▶ Stays put while remaining easy to remove



EndoSequence® Flowable Dam Material

4	3mL Syringes
40	Tips
1	IFU

Retail: \$104.95

Special: \$73.47

ENDOSEQUENCE® BC LINER™

- ▶ RMGI optimized for use with BC Sealer/RRM
- ▶ Self-bonding and dual-cure
- ▶ Strong (241 mPA) and versatile
- ▶ Excellent for orifice barriers



Now
available
in blue!



EndoSequence® BC Liner™ (White or Blue)


1	BC Liner
40	BC Liner Short Tips
1	Safety Data Sheet (SDS)

Retail: \$129.95

Special: \$99.95


ESR™ CM RECIPROCATING FILES

ESR CM utilizes the latest NiTi technology to provide clinicians with superior safety, efficiency and performance. ESR CM files are also designed to preserve more coronal tooth structure compared to WaveOne® Gold.




HIGHER FLEXIBILITY

UP TO **X3.3***




GREATER RESISTANCE TO CYCLIC FATIGUE

UP TO **X3.6***




LOWER SCREWING EFFECT

UP TO **42%***



LOWER STRESS ON DENTIN

UP TO **40%***



**Text results available upon request*

ESR CM files are provided in sterile blister packs for optimal convenience and safety. If sterility is compromised prior to treatment the files can be sterilized without any loss of functionality or performance.



File/Tip ID	Tip Size	21 mm	25 mm	31 mm
△ White	15	5028076U0	5028084U0	5028092U0
▲ Red	25	5028078U0	5028086U0	5028094U0
▲ Blue	30	5028079U0	5028087U0	5028095U0
▲ Black	40	5028081U0	5028089U0	5028097U0
▲ Yellow	50	5028082U0	5028090U0	5028098U0


ESR™ CM Files (3/Pack)

Retail: \$61.50/Pk

Buy 50 Packs, Get 40% Off!
Buy 25 Packs, Get 25% Off!
Buy 10 Packs, Get 15% Off!


Mix and match

RECIPROICATION REDEFINED...



- Reverse-cutting flute design: cuts on the counter-clockwise motion
- Compatible with EndoSync™ Plus, ProMark®, e3® and NLZ™ E electric motors
- Available in #15 (Scout), #25, #30, #40 and #50
- Heat treated NiTi allows for pre-bending
- Constant taper design
- Sterile blister packs of three instruments

STERILE R




ESR™ CM VS. WAVEONE® GOLD

- ESR CM features a similar but more flexible flute design*

**63% of evaluators noted that ESR CM is more flexible than WaveOne Gold in an internal clinical survey.*


- ESR CM features a size #30 which is less invasive than the WaveOne® Gold primary file
- ESR CM can be autoclaved if sterility is compromised
- ESR CM features a novel safety tip optimized for reciprocation

WaveOne® Gold #25 Primary



Modified Rectangular Cross-Section


ESR CM #30



Modified Triangular Cross-Section


NLZ™ is a trademark of Nakanishi, Inc.
WaveOne®, ProMark® and e3® are registered trademarks of Dentsply Sirona, Inc.

ENDOSEQUENCE™ CM TAPER CONTROLLED MEMORY FILES



Uncompromised Quality. Now in a Progressive Taper.

Introducing the latest addition to our popular EndoSequence CM offering. The EndoSequence CM taper file provides clinicians with a progressive tapered instrument that delivers the cutting efficiency and flexibility of EndoSequence CM in the familiar shapes of ProTaper Gold™ files.



File/ Tip ID	Tip Size	Variable Taper		
		16 mm	21 mm	25 mm
▲ N/A	SX/19	5027891U0	—	—
▲ Purple	S1/18	—	5027892U0	5027899U0
△ White	S2/20	—	5027893U0	5027900U0
▲ Yellow	F1/20	—	5027894U0	5027901U0
▲ Red	F2/25	—	5027895U0	5027902U0
▲ Blue	F3/30	—	5027896U0	5027903U0
Assortment	SX-F3	—	5027898U0	5027905U0
▲ Black	F4/40	—	5027897U0	5027904U0

EndoSequence® CM Files (4/Pack)

Retail: \$60.90/Pk

Buy 50 Packs, Get 40% Off!

Buy 25 Packs, Get 25% Off!

Buy 10 Packs, Get 15% Off!

Mix and match

ProTaper Gold™ is a trademark of Dentsply Sirona, Inc.

ENDOSEQUENCE® BC TEMP™

BC Temp is a premixed bioceramic paste for intracanal dressing. Unlike traditional Calcium Hydroxide pastes, the primary components of BC Temp are Calcium Silicates and Calcium Oxide which combine with the water naturally present in dentin to produce Calcium Hydroxide which dissociates into Ca^{2+} and OH^- . The hydroxyl ions (OH^-) released are responsible for a significant increase in the pH of the surrounding tissues, making the environment unsuitable for bacterial growth. The hydration reaction of Calcium Silicates allows BC Temp to have a gradual and slow release of calcium and hydroxyl ions which eliminates the need for frequent applications.

BC Temp is highly radiopaque (>9mm) and utilizes Calcium Tungstate as a radiopacifier so it is non staining and unlike conventional Calcium Hydroxide materials, BC Temp is easily removed which reduces chair time.

- ▶ Single or Less Frequent Application
- ▶ Highly Biocompatible
- ▶ Optimal pH
- ▶ Easily Removed
- ▶ High Radiopacity
- ▶ Non-staining
- ▶ Optimal flow – $23,85 \pm 0,33$ mm




EndoSequence® BC Temp™

4 .5g Syringes
20 Tips
1 IFU


Retail: \$70.00

Special: \$50.50



Initial X-Ray

Radiograph after initial placement of BC Temp



Final X-Ray

Final obturation; 15 days after initial placement of BC Temp

*Images courtesy of Dr. Kely F. Bruno
Universidade Federal de Goiás, Goiânia, GO, Brazil*

Test results available upon request

BRASSELER U.S.A. DENTAL, LLC

One Brasseler Blvd.
Savannah, GA 31419

Toll Free: 800.841.4522

Phone: 912.925.8525

Fax: 888.610.1937

Shop.BrasselerUSA.com



BRASSELER
USA® BY YOUR SIDE

©2021 Brasseler USA. All rights reserved. Printed in the U.S.A.

B-5500-US-55M-07.21-RRD



BRASSELER
USA® RESTORATIVES

***Innovation Meets
Excellence...***

*A Partnership for
Restorative Success*

Coming in August 2021!



powered by
**Clinician's
Choice**

Invoice or statement prices may reflect or be subjected to a bundled discount or rebate pursuant to purchase offer, promotion, or discount program. You must fully and accurately report to Medicare, Medicaid, Tricare and/or any other federal or State program, upon request by such program, the discounted price(s) or net price(s) for each invoiced item, after giving effect to any applicable discounts or rebates, which price(s) may differ from the extended prices set forth on your invoice. Accordingly, you should retain your invoice and all relevant information for your records. It is your responsibility to review any agreements or other documents, including offers or promotions, applicable to the invoiced products/prices to determine if your purchase(s) are subject to a bundled discount or rebate. Any such discounts must be calculated pursuant to the terms of the applicable purchase offer, promotion, or discount program. Participation in a promotional discount program is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the discount program requirements.

Distributed by Brasseler USA