


stability where it counts



Tapered Immediate Molar

Tapered IM implants offer a solution for molar extraction sockets. Available in 7 and 8mm diameters, the implant design features deep aggressive threads for compressive bone loading and primary stability. The reduced collar diameter ensures the implant does not deviate from its intended placement location in extraction sites. The platform-switched, dual-affinity Laser-Lok® surface offers crestal bone retention and connective tissue attachment for flexible placement.

For more information, contact BioHorizons
Customer Care: (phone number) or
(shop or visit us online) at www.biohorizons.com

BIOHORIZONS®
SCIENCE • INNOVATION • SERVICE