

MedLED®
AAOMS Meeting Sale
Valid September 30 - October 2, 2021



MedLED MC7
\$3,995.00 US
\$4,594.25 CA

Purchase the **MedLED MC7** Standard Kit,
receive a free upgrade to the Hospital Kit.

S

Standard Kit

ML-MC7-HT-SK
ML-MC7-ST-SK



- 1 MedLED Headlight
- 2 LiPO Batteries
- 1 USB Wall Charger Set
- 1 Carrying Case & Bag

+

Hospital Kit

add ML-M-90611-P



- 1 MedLED Headlight
- 6 LiPO Batteries
- 1 Battery Charger, 4-Bay
- 1 USB Wall Charger Set
- 1 Carrying Case & Bag

To purchase, email contact@klsmartin.com

KLS martin®
GROUP

Products shown may not be licensed or cleared in all markets.
Availability is subject to change without notice. Licensed for US and CA.

MedLED®
AAOMS Meeting Sale
 Valid September 30 - October 2, 2021



MedLED MC6
 \$2,995.00 US
 \$3,444.25 CA

MedLED MC3
 \$1,750.00 US
 \$2,012.50 CA

Purchase the **MedLED MC3** and / or **MC6** Standard Kit,
 receive a free upgrade to the Hospital Kit.

S

Standard Kit

ML-MC3-HT-SK, ML-MC3-ST-SK
 ML-MC6-HT-SK, ML-MC6-ST-SK



1 MedLED Headlight
 2 LiPO Batteries
 1 USB Wall Charger Set
 1 Carrying Case & Bag

+

Hospital Kit

for MC3, add ML-M-90612-P
 for MC6, add ML-M-90611-P



1 MedLED Headlight
 6 LiPO Batteries
 1 Battery Charger, 4-Bay
 1 USB Wall Charger Set
 1 2x Snap-on Magnifiers
 1 Carrying Case & Bag

To purchase, email contact@klsmartin.com

KLS martin®
 GROUP

KLS-Martin L.P.
 P.O. Box 16369 · Jacksonville, FL 32245-6369
 11201 Saint Johns Industrial Pkwy S · Jacksonville, FL 32246-7652
 904.641.7746 · 800.625.1557 · Fax 904.641.7378
usa@klsmartin.com · www.klsmartinnorthamerica.com

91-600-06-07 · REV 06 · 2021.07 · medLED AAOMS Promotion · Offer good for a limited time only on retail prices. Limit 10 per customer. Discount per purchaser or affiliated purchasers. Not valid with any other offers or discounts. Products shown may not be licensed or cleared in all markets. Availability is subject to change without notice. The information presented is intended to demonstrate a KLS-Martin L.P. product, please refer to the IFU for the full list of indications, contraindications, warnings, precautions, and sterilization information. Always refer to the IFU before using any KLS-Martin L.P. product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients.

The top half of the image features a dark blue background. In the upper left, the KLSmartin Group logo is displayed, with 'KLS' in red and 'martin' in white, followed by 'GROUP' in smaller white capital letters. To the right of the logo, a medical headlight device is shown. It is a black, cylindrical unit mounted on a bracket, which is attached to a larger black structure with a grid of small holes. The headlight is angled downwards and to the right, emitting a bright, focused beam of light that illuminates the dark blue background.


KLSmartin[®]
GROUP

MedLEDChrome[®]

Truly Cordless[™] Headlight Technology

www.klsmartinnorthamerica.com

Products shown may not be licensed or cleared in all markets.
Availability is subject to change without notice. Licensed for US and CA.



Concentrate on the brightly lit surgical field with total freedom and complete comfort.

MedLED Chrome® Headlights

Freedom without sacrifice

The MedLED Chrome® headlight has incorporated the best that technology offers to accomplish the once impossible: a silent, comfortable, and cordless headlight that can light up any field in the O.R. in terms of both intensity and spot size.

Features

- Luminous intensity:
MC3 – 100,000 lux
MC6 – 200,000 lux
MC7 PRO – 275,000 lux
- Battery service life (infinite with 4-bay charger):
MC3 – 7.75 hr
MC6 – 6.75 hr
MC7 PRO – 3.50 hr
- USB charging time: 6.5 hours
- 4-Bay charging time: 10 hours; charge up to two sets at a time
- Adjustability: full beam angle adjustability, left-to-right swivel, increased spot size range
- Spot quality: crisp, circular beam of homogenous light


MedLED Chrome®

Freedom without sacrifice

If you have a need for a powerful, long lasting, comfortable headlight with no cord, KLS Martin has you covered with a full range of headlights and accessories.



Power and Brightness



Spot Adjustability and Control

200 mm
MC 7 PRO

Maximum spot diameter range
1:1 scale

100 mm
MC 3


90 mm
MC 6

MedLED Chrome®

Head-Mounted Circuit Technology

Lose the optical cables and belt pack cords and go cordless without sacrificing basic operating light needs.

Does not interfere with sterility. Without cables or battery belt packs, it's kept safely above the neck.




Quiet-Cool™ Heatsink Technology

Cool and quiet


Proprietary technology is used to cool the LED instantly, avoiding the potential for unnecessary heat buildup. We also designed the housing itself to manage any remaining thermal dissipation.

Fans demand more power, requiring large batteries attached to the waist of the user. As a moving part in constant operation, the fan motor has a high potential for malfunction.


MC 3






MC 6



MC 7 PRO



MedLED Chrome® Headlights Ready to Use!


	<p>MedLED Chrome® MC 3</p> <p>Consisting of:</p> <ul style="list-style-type: none"> • Chrome® MC 3 headlight • Carrying case with shoulder strap • Universal wall charger • 2 Rechargeable lithium-polymer batteries 	<p>Item No.</p> <p>ML-MC3-HT-SK ML-MC3-ST-SK ML-MC3-HT-HK* ML-MC3-ST-HK*</p>
	<p>MedLED Chrome® MC 6</p> <p>Consisting of:</p> <ul style="list-style-type: none"> • Chrome® MC 6 headlight • Carrying case with shoulder strap • Universal wall charger • 2 Rechargeable lithium-polymer batteries 	<p>Item No.</p> <p>ML-MC6-HT-SK ML-MC6-ST-SK ML-MC6-HT-HK* ML-MC6-ST-HK*</p>
	<p>MedLED Chrome® MC 7 PRO</p> <p>Consisting of:</p> <ul style="list-style-type: none"> • Chrome® MC 7 PRO headlight • Carrying case with shoulder strap • Universal wall charger • 2 Rechargeable lithium-polymer batteries 	<p>Item No.</p> <p>ML-MC7-HT-SK ML-MC7-ST-SK ML-MC7-HT-HK* ML-MC7-ST-HK*</p>

ST stands for soft top, HT stands for hard top


* Hospital Kit includes 4 bay charger and additional batteries.

MedLED Chrome® Accessories


Magnifier lenses MC 3

	Diopter	Magnification	Increase in eyesight	Optical power	Recommended working distance (focal length)	Item No.
	+ 2	+ 0.40	40 %	1.4 x	250 mm	ML-M-90301
	+ 3	+ 0.60	60 %	1.6 x	220 mm	ML-M-90303
	+ 4	+ 0.80	80 %	1.8 x	200 mm	ML-M-90305

Magnifier lenses MC 6 / MC 7 PRO


	Diopter	Magnification	Increase in eyesight	Optical power	Recommended working distance (focal length)	Item No.
	+ 2	+ 0.40	40 %	1.4 x	250 mm	ML-M-90321
	+ 3	+ 0.60	60 %	1.6 x	220 mm	ML-M-90323
	+ 4	+ 0.80	80 %	1.8 x	200 mm	ML-M-90325

Desktop charger



ML-M-90202


Spare batteries, rechargeable (2 pcs.)




ML-M-90004

MedLED Chrome® Accessories

Spare head bands




Spare headstrap




Bezel




GoPro Adapter




Amber filter



Desktop charger/Mains plug



*Refer to the Instructions For Use for proper cleaning and maintenance.



truly cordless™

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

**KLS Martin Medical (Shanghai)
International Trading Co. Ltd.**

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Tuttlingen · Germany
Tel. +49 7461 706-0
info@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 505 7838
malaysia@klsmartin.com

KLS Martin de México S.A. de C.V.

Mexico City · Mexico
mexico@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

London · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS-Martin L.P.**A company of the KLS Martin Group**

Mail: P.O. Box 16369, Jacksonville, FL 32245-6369
Office: 11201 St. Johns Ind Pkwy S., Jacksonville, FL 32246-7652
Tel. 904.641.7746 · 800.625.1557
Fax 904.641.7378
usa@klsmartin.com · www.klsmartinnorthamerica.com

ORAL SURGERY
motor system

IMPLANTOLOGY
motor system

BienAir⁺
Dental

DISTRIBUTED BY

KLS martin® L.P.

NEW OMS

- INTEGRAL COOLING FAN
- HIGH TORQUE
- SPEED STABILITY



TRY OMS WITH 1:2 STRAIGHT HANDPIECE,
DESIGNED SPECIFICALLY FOR IMPACTED THIRD MOLAR.



The high rotational speed reduces cutting time
by **70%** (12 seconds to completely section a tooth),
while significantly reducing the force required.

12 sec

80,000 rpm

40 sec

40,000 rpm

S W I S S  M A D E

MICROMOTOR MX-i LED PLUS



HIGH TORQUE

The system allows up to 80 Ncm torque with 20:1 CA.

SPEED STABILITY

High Torque quantity does not equal stable mechanical power. **OMS** is the most stable motor in low RPM (implant placement) ever measured, providing the user with confidence unmatched by competition.

INTEGRAL COOLING FAN

MX-i LED PLUS micromotor supports use of HP 1:2 for high-speed impacted 3rd Molar (Wisdom Tooth) extraction.



ADVANTAGES OF OMS SYSTEM

Recommended procedures:

- Oral Surgery (including impacted third molar extraction) ✓
- Implant placement ✓

Max. Torque in implant mode	80 Ncm
-----------------------------	--------

Torque control	✓
----------------	---

Thermo-washer disinfectable (MX-i LED PLUS and cable)	✓
---	---

Integral cooling fan (MX-i LED PLUS)	✓
--------------------------------------	---

LED Light (MX-i LED PLUS)	✓
---------------------------	---

Supported handpieces and contrangles:

- Contra-angle: CA 20:1, CA 1:2.5 ✓
- Handpieces: PM 1:1, PM 1:2 ✓

S W I S S  M A D E

FOR MORE INFORMATION:

Contact your local KLS-Martin Sales Consultant
or call the corporate office at 800.625.1557

KLS-Martin L.P.

A company of the KLS Martin Group

Mail: P.O. Box 16369, Jacksonville, FL 32245-6369

Office: 11201 St. Johns Ind Pkwy S., Jacksonville, FL 32246-7652

Tel. 904.641.7746 · 800.625.1557

Fax 904.641.7378

usa@klsmartin.com · www.klsmartinnorthamerica.com