




Your in-lab
ceramic solution

neoss® | NeoBase™ Abutments

Intelligent Simplicity

A complete range of titanium base abutments for local milling of ceramic restorations with angulated and straight screw channels.



Full flexibility for individual needs

- Variety of post heights for additional retention area
- Variety of gingival margin heights for ideal implant placement
- All NeoBase™ abutments cater for both straight and angulated screw channels

Angulated screw channels – improved handling and esthetics

- Favorable screw channel exit positions
- Screw-retained solution in the anterior for improved esthetics and occlusal function
- Easier posterior access

Design freedom and production adaptability

- For use at angles up to 25° with 360° rotational freedom
- Compact inner and outer geometry for maximum design freedom and compatibility with conventional milling machines
- NeoBase™ post is designed for reliable retention and maximum ceramic thickness



Neoss provides a complete range of products to support your digital workflow and enable scanning of implant positions, design of prosthetics, and alternatives for local and central production.

3shape  exocad

Design software and Libraries

Neoss provides design libraries for 3Shape and Exocad software to enable you to create your digital prosthetic design on NeoBase™ abutments.

Libraries together with the installation guidelines are available for download on www.neoss-cadcam.de.

Design

A range of design options are available, including both single and multiple ceramic restorations, as well as the Esthetiline solution with matching gingival design and anatomical fit from healing abutments to final restorations.

Milling

NeoBase™ mating geometry, including retentive features, is designed to maximize the options of compatible dental ceramic milling machines.

Final result

NeoBase™ is the foundation for esthetically driven retrievable restorations, avoiding the need for intra-oral cementation. NeoBase™ is your local solution for ceramic restorations without compromise on design.



Neoss has been providing stress-free ceramic restorations for over 15 years with the original NeoLink® concept. The new range of NeoBase™ abutments caters for both straight and angulated screw channels, and is yet another example of our unique **Intelligent Simplicity** in action. Innovative and cost-effective solutions putting you, your team and your patients first.

Assortment

Art. No. Description

31367	NeoBase™ Mono G0.3 mm – H3.6 mm
31368	NeoBase™ Mono G1.5 mm – H3.6 mm
31369	NeoBase™ Mono G0.3 mm – H5.6 mm
31370	NeoBase™ Mono G1.5 mm – H5.6 mm
31371	NeoBase™ Multi G0.3 mm – H3.6 mm
31372	NeoBase™ Mono G0.3 mm – H3.6 mm Ø3.25
31373	NeoBase™ Mono G0.3 mm – H5.6 mm Ø3.25
31374	NeoBase™ Multi G0.3 mm – H3.6 mm Ø3.25

Contact your local representative for information about availability in your market.

Contact details

Australia

Neoss Australia Pty. Ltd
PO Box 404
New Farm
QLD 4005
T +61 7 3216 0165
F +61 7 3216 0135
E info.au@neoss.com

Austria

Neoss GmbH
T +49 221 55405-322
F +49 221 55405-522
E info@neoss.de

China

Shanghai Legend Life Technology Co., Ltd
4th floor, Block 6
1525 Minqiang Road
Songjiang District
Shanghai 201612
P. R. China
T +86-21-6765 7130-617
F +86-21-6765 7131
E info@legendlifetech.com

Croatia

Medical Intertrade d.o.o
10431 Sveta Nedelja
HRVATSKA
T +385 1 33 74 020
F +385 1 33 36 035
E dental@medical-intertrade.hr

Denmark

Neoss AB
T +46 (0)31 88 12 80
E info@neoss.se

Germany

Neoss GmbH
Im MediaPark 8
DE-50670 Köln
T +49 221 55405-322
F +49 221 55405-522
E info@neoss.de

Italy

Neoss Italia S.r.l.
Viale Certosa 138
IT-20156 Milano
T +39 02 92952 1
F +39 02 92952 250
E italia@neoss.info

Japan

Neoss Japan Co.
Y101 Building 5F
No. 1-4 Kanda-Sudacho Chiyodaku
Tokyo 101-0041
T +81 3 5289 3511
F +81 3 5289 3512
E info@neoss.co.jp

Middle East & Africa

Neoss Limited
Unit 3737 Level 1
Jewelry & Gemplex 3
P.O. Box 5003315
Dubai
United Arab Emirates
E info.me@neoss.com

Netherlands

Rotor, Medical Supplies
Economiestraat 16
NL-6135 KV Sittard
T +31-46-481 12 03
F +31-46-481 12 04
E mail@rotor-medical.com

New Zealand

Neoss Australia Pty. Ltd
T +61 7 3216 0165
F +61 7 3216 0135
E info.au@neoss.com

Norway

Neoss AB
T +46 (0)31 88 12 80
E info@neoss.se

Republic of Ireland

Neoss Ltd
T +44 1423 817-733
F +44 1423 817-744
E info@neoss.com

Spain

Formación, Gestión Y Rendimiento, S.L
C. Londres 54 E2B
ES-08036 Barcelona
T +34 931 408 770
E info@fgrdental.com

Sweden

Neoss AB
Arvid Wallgrens backe 20
SE-413 46 Göteborg
T +46 (0)31 88 12 80
E info@neoss.se

Arc Solutions AB

Box 13011
SE-250 13 Helsingborg
T +46 (0)42 3017440
E info@arcsolutions.se

Switzerland

Flexident AG
Amlehnstr. 22
CH-6010 Kriens
T +41413104020
F +41413104025
E info@flexident.ch

Thailand

S.D.Tantavade Co., Ltd.
47 Samakkee 20, Samakkee Road
Thasai Muang
Nonthanburi 11000
T +66 2 952 4815 18
F +66 2 952 4819
E info@sdt1988.com

Turkey

Dentmark Sağlık Ürünleri Anonim Şirketi
Girne Mah. Girne Cad. No: 117
Omega Deniz Plaza Kat: 2/12
Maltepe
İstanbul
T +90 216 417 4040
F +90 216 417 4041
E info@neossshi.com

United Kingdom

Neoss Ltd
Windsor House
Cornwall Road
Harrogate
HG1 2PW
T +44 1423 817-733
F +44 1423 817-744
E info@neoss.com

United States

Neoss Inc.
21860 Burbank Blvd.,
Ste. 190 (North Entrance)
Woodland Hills, CA 91367
T +1 866 626-3677
F +1 818 432-2640
E contact.usa@neoss.com