

ReStore[™]

The ReStore Soft Exo-Suit A Revolution in Post-Stroke Gait Training

Functional

The ReStore soft design combines natural movements with plantar-flexion and dorsiflexion assistance that adaptively synchronize with the patient's own gait to facilitate functional gait training.

Versatile

Individualized assistance and compatibility with supplemental support aids ensure that ReStore has broad applications for patients across the gait rehabilitation spectrum.

Data-Driven


Real time feedback and adjustable levels of assistance enable the therapist to optimize sessions and track each patient's progress.





ReWalk[®]

ReStore

How Does ReStore Work?



Sensors on each foot detect gait events and synchronize timing of assistance



Posterior cable contracts to provide PF assistance for **forward propulsion**

Anterior cable contracts to provide DF assistance for **ground clearance**

Slack and Brace modes enable additional gait training capabilities

ReStore App provides real time feedback and session control to optimize results



Seamless integration with a wide array of functional gait training activities within your clinic

To learn more, visit us at www.rewalk.com or contact us at www.rewalk.com/contact
ReWalk Robotics • 200 Donald Lynch Boulevard, Marlborough, MA 01752



© 2019 ReWalk Robotics, Inc. MK40R01 Rev.02