

# LAPspay® Systems for Dogs and Cats

## Laparoscopic Ovariectomy / Ovariohysterectomy

### Single and Multi-Port



# LAPspay® Multi-Port

## Standard Laparoscope



- 62046AA **HOPKINS® Straight Forward Telescope 0°**, enlarged view, diameter 5 mm, length 29 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green

## TERNAMIAN ENDOTIP Cannulas



- 60120MS TERNAMIAN **ENDOTIP Cannula**, with thread and without insufflation stopcock, size 6 mm, working length 6.5 cm, color code: black-white including: **Silicone Leaflet Valve**

- 60120MSS **Same**, with insufflation stopcock


## Standard Trocar and Cannula, with multifunctional valve



- 62160WZ **Trocar**, size 6 mm, color code: black-white including: **Trocar**, only, with pyramidal tip **Cannula**, without valve, with insufflation stopcock, working length 8.5 cm **Multifunctional Valve**


## LAPspay® Single Port

### Operating Laparoscope



- 26034AAV **HOPKINS® Operating Telescope 0°**, with parallel eyepiece, diameter 10 mm, length 23 cm, **autoclavable**, fiber optic light transmission incorporated, with 6 mm working channel, including adaptor with Silicone Leaflet Valve 26034V

### TERNAMIAN ENDOTIP Cannula, with multifunctional valve



- 60103TMR TERNAMIAN **ENDOTIP** Cannula, with thread and rotatable insufflation stopcock, size 11 mm, working length 10.5 cm, color code: green  
including:  
**Multifunctional Valve**  
**Seal Set 60/10**  
**Mandrel**, for opening automatic valves

### Standard Trocar and Cannula



- 62103FZ **Trocar**, size 11 mm, color code: green  
including:  
**Trocar**, only, with pyramidal tip  
**Cannula**, without valve, with insufflation stopcock, working length 8.5 cm  
**Automatic Valve**



- 62141DB **Reducer**, 11/5 mm

**Whether you prefer the Multi-Port or Single-Port technique, minimally invasive ovariectomy and ovariohysterectomy using a KARL STORZ LAPspay® surgery system offers the following advantages:**

- Unsurpassed optics for a clear, bright, enlarged view
- Less pain for the patient since tissues are cauterized and transected with precision, rather than torn
- Complications, such as wound dehiscence, ovarian remnants and ligated ureters, are reduced to a minimum
- Reduced chance of infection due to minimal exposure of the abdomen
- Faster recovery due to small incision size
- Satisfied clients who demand the best for their pets
- Practice expansion through state-of-the-art surgical services
- The opportunity to improve your laparoscopy skills, which can be applied to numerous other minimal access surgical techniques



Grasping Forceps, 5 mm



66322MG

**CLICKLINE Grasping Forceps**, “tiger jaws”, rotating, dismantling, insulated, with connector pin for unipolar coagulation, with LUER-Lock irrigation connector for cleaning, single action jaws, 2 x 4 teeth, for grasping and removal of solid organs, size 5 mm, length 36 cm including:

**Plastic Handle**, with style ratchet

**Metal Outer Sheath**, insulated

**Forceps Insert**

33422MG

**Same**, length 43 cm

Bipolar Forceps, 5 mm



66351DN

**ROBI PLUS®**, rotating, with connector pin for bipolar coagulation, with integrated cutting function, single action jaws, size 5 mm, length 35 cm including:

**Ring Handle**, bipolar, with activator for cutting function

**Outer Sheath**, insulated, length 35 cm

2x **Knife Sheath with Integrated Blade**, recommended for single use, length 35 cm

**Forceps Insert**, length 35 cm

66451DN

**Same**, length 42 cm

Single use knife sheaths for ROBI PLUS®



66300M

**Knife Sheath with Integrated Blade**, recommended for single use, length 35 cm, for use with Bipolar Instrument ROBI PLUS® 66351DN

66400M

**Same**, length 42 cm, for use with Bipolar Instrument ROBI PLUS® 66451DN

## Additional Instruments



- 62175H TANKERSLEY **Ovariectomy Hook**, with heavy handle, for percutaneous use, working length 8.5 cm
- 62175HK **Same**, with small hook, for small dogs and cats



- 62120J VERESS **Pneumoperitoneum Needle**, with spring-loaded blunt stylet, LUER-Lock, diameter 2.1 mm, length 10 cm

## Insufflator



- 62432501 **CO<sub>2</sub> ENDO-ARTHROFLATOR-VET**, power supply 100-240 VAC, 50/60 Hz including:  
**Silicon Tubing Set**, sterilizable  
**Universal Wrench**  
**Gas Filter**, sterile, for single use, package of 10

**The pressure hose and CO<sub>2</sub> bottles must be ordered separately, see catalog *Veterinary Endoscopy, Small Animals***

## Electrosurgery

### AUTOCON® III 300

The AUTOCON® III 300 is a versatile electrosurgical unit that can be used in bipolar or unipolar mode for a variety of surgical procedures. It delivers precise coagulation and cutting in both endoscopic and open procedures.

KARL STORZ offers an assortment of instruments for use with the AUTOCON® III 300, for use in a wide range of veterinary applications.



UH300

**AUTOCON® III 300**, power supply 200-240 VAC, 50/60 Hz, including mains cord

**HF connecting sockets unipolar:**

1x 3-pin US type

1x 5 mm connector KARL STORZ/Erbe VIO

**HF connecting sockets bipolar:**

1x 2-pin US type (28.58 mm)

2x KARL STORZ/Erbe VIO

**Connecting socket:**

Neutral electrode 2-pin

Can be integrated into OR1™ with KARL STORZ OR1™ SCB CONTROL WU 300

System requirements: SCB R-UI Software Release 20090001-46 or higher

UH300U

**Same**, power supply 100-127 VAC

### High Frequency Cord for use with ROBI PLUS® and AUTOCON® III 300



26176LE

**Bipolar High Frequency Cord**, length 300 cm, for AUTOCON® III 300/400, AUTOCON® II 400 SCB system (111, 113, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, Coagulator 26021B/C/D, 860021B/C/D, 27810B/C/D, 28810B/C/D, AUTOCON® series (50, 200, 350), Erbe-Coagulator, T and ICC series



UF902

**Two-Pedal HF Footswitch**, with button for switchover function, for use with HF generators

See catalog, *Veterinary Endoscopy, Small Animals*, for a full range of instruments and high frequency cords for use with the AUTOCON® III 300



Register now to start receiving FREE client referrals for endoscopy procedures.



# 75 Years

*Shaping the Future  
of Endoscopy with you*

**STORZ**  
KARL STORZ — ENDOSKOPE  
THE DIAMOND STANDARD

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [vet@karlstorz.com](mailto:vet@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)

