

# Cusp-Lok

## Cuspid Traction System

A surgically placed orthodontic appliance, The Cusp-Lok family of brackets solves the problem of impacted cuspids for both the surgeon and orthodontist. Each bracket is manufactured with 14k gold chain, either 3/4" or 1" in length, with optimal aperture for ligature wire or power thread.

**Note:** Quantity discounts available: 10% for 3-4 boxes, and 15% 5+ boxes

### Regular/Low Profile

The original and still most preferred solid bracket. Offers a medium curved pad, and is best for applications with the facial surface of cuspids. As with all brackets, the 14k chain is hypo-allergenic and kind to tissue.



CLI &  
CLII

Code#	Profile/Description	Chain Length	5pk	10pk
CLI	CUSP-LOK I	3/4"	\$140	\$240
CLII	CUSP-LOK II	1"	\$150	\$255
CLILP	CUSP-LOK I LP	3/4"	\$140	\$240
CLILP	CUSP-LOK II LP	1"	\$150	\$255



CLILP &  
CLILP

### Mini Flat/Low Profile

Ideal for lingual placement, or teeth other than cuspids. As with other styles, the low profile brackets have the first link of chain attached parallel to the bonding pad, eliminating the "bump" of tissue in these cases.



CLIM &  
CLIIM

Code#	Profile/Description	Chain Length	5pk	10pk
CLIM	CUSP-LOK I Mini	3/4"	\$155	\$255
CLIIM	CUSP-LOK II Mini	1"	\$170	\$275
CLIMLP	CUSP-LOK I Mini LP	3/4"	\$155	\$255
CLIIMLP	CUSP-LOK II Mini LP	1"	\$170	\$275



CLIMLP &  
CLIIMLP

### Mini Curved/Low Profile

Ideal for cases with deep impaction where limited facial exposure can be expected. Similar in shape to the Cusp-Lok I brackets, these brackets have 25% less surface area and rounded corners for placement in difficult to reach areas.




CLIMC &  
CLIIMC

Code#	Profile/Description	Chain Length	5pk	10pk
CLIMC	CUSP-LOK I Mini	3/4"	\$155	\$255
CLIIMC	CUSP-LOK II Mini	1"	\$170	\$275
CLIMCLP	CUSP-LOK I Mini LP	3/4"	\$155	\$255
CLIIMCLP	CUSP-LOK II Mini LP	1"	\$170	\$275



CLIMCLP &  
CLIIMCLP

### Cusp-Lok Placement Instruments



## Chain & Mesh/Low Profile

Designed to permit bonding agents to flow through the 4mm diameter stainless mesh, resulting in improved bonding strength. Mesh pad can be manually contoured to any shape quickly and easily. Works very well with our Light Cure bonding agent.



CLCMR

Code#	Profile/Description	Chain Length	5pk	10pk
CLCMR	CUSP-LOK Chain & Mesh	1"	\$150	\$255
CLCMLP	CUSP-LOK Chain & Mesh LP	1"	\$150	\$255



CLCMLP

## Swivel 360

Swivel bracket has no solder joint, permitting activation in any direction for direct vector force. Lingual button pad measures 3mm in diameter.



CL360

Code#	Profile/Description	Chain Length	5pk	10pk
CL360	CUSP-LOK Swivel 360	1"	\$145	\$235

## Molar Bracket

Features a molar sized, curved pad with contours. For maximum surface area and bonding strength.



CLMB

Code#	Profile/Description	Chain Length	5pk	10pk
CLMB	Molar Bracket	1"	\$165	\$270
CLMBLP	Molar Bracket LP	1"	\$165	\$270



CLMBLP

## Titanium Spring

Suitable for vertically impacted teeth when long linear movements are required. Spring is crafted with nickel titanium.



CLTS

Code#	Profile/Description	Length	Cost Per Piece
CLTS	Titanium Spring Bracket	1/2" Spring	\$40

## Cusp-Lok Complete Kit

The Cusp-Lok Complete Kit contains everything necessary for an expose and bond procedure.



### CLKTSC

Cusp-Lok Complete Kit  
(custom kits upon request)  
**\$345**

### Includes:

- (2) Cusp-Lok I
- (2) Cusp-Lok I Mini
- (2) Cusp-Lok I Low Profile
- (2) Cusp-Lok II
- (2) Cusp-Lok Chain & Mesh
- (2) Cusp-Lok Chain & Mesh LP
- (1) Placement Instrument
- (1) Crown Lok A & B Paste
- Etchant Gel and Accessories

## Quick, Effective Bonding Without a Curing Light

Self Cure bonding agent is the perfect system for those who desire a quick and uncomplicated bonding solution, no curing light required. Achieves a solid bond 2 minutes after application.

### CLSC

Self Cure bonding kit  
**\$70**

## Light Cure Package

- Cordless Curing Light
- Light Cure Kit
- 10 Chain & Mesh Brackets

### LCP

Light Cure Package  
**\$550**



## Rapid Light Curing Results

Our Light Cure kit is a single paste light cure bonding agent that is designed for a wet or dry environment. When used with Cusp-Lok Chain & Mesh brackets, it permits the bonding agent to be completely cured, instead of just along the perimeter of the bracket.

### CLLC

Light Cure bonding kit  
**\$115**

[www.xemax.com](http://www.xemax.com)

**XEMAX**  
SURGICAL PRODUCTS, INC.

**1-800-257-9470**