



X-OTOMES

by A. Titan

ATRAUMATIC EXTRACTIONS

An Instructional Tip: Beveled tips follow root structure down to PDL. Using a 2-3 degree rotation around circumference of tooth you will gently separate the fibers from the root structure down to the PDL, creating a trough that will minimize tissue and bone trauma.



XO-10

Small straight tip works on anterior teeth



XO-20

Small curved tip works on anterior and single rooted teeth



XO-30

Small angled tips work on the mesial of anteriors & bicusps



XO-40

Small angled tips work on the distal of anteriors & bicusps



XO-50

Medium angled tips work on 1st, 2nd, and Premolars



XO-60

Medium angled tips work on 1st, 2nd, and Premolars



XO-80

Straight Micro



XO-110

Notched tip X-Otome



XO-70

Straight spade shaped blade ideal for interproximal tissue reflecting



XO-77

Beveled spade shaped blade slides nicely below gum line and follows the contour of root to both sever fibers and aid in root removal



XO-77D

Angled for the distal, the beveled spade shaped blade slides nicely below gum line and follows the contour of root to both sever fibers and chase ankylosed roots



XO-77M

Angled for the mesial, the beveled spade shaped blade slides nicely below gum line and follows the contour of root to both sever fibers and chase ankylosed roots



XO-77AL

Angled left, the beveled spade shaped blade slides below gum line to sever fibers and chase ankylosed roots



XO-77AR

Angled right, the beveled spade shaped blade slides below gum line to sever fibers and chase ankylosed roots



XO-90

Large angled tip excellent for 3rd molars




XO-100

Large angled tip gives the ability to maneuver and elevate tuberosity around 3rd molar

- Made in the America
- Hybrid instrumentation, a cross between an Elevator and a Periotome
- Custom sure-grip handle for controlled luxation force
- Honed narrow tips to disrupt the periodontal fibers
- Preserves tissue structures
- Eliminates bone fracture and promotes minimally invasive exodontia
- 100% stainless steel

Elemental Kit K-3300

Contains:
XO-10, XO-20, XO-77 and
a sterilization cassette




X-Otome Kit K-6600

Includes one of each:
XO-10, XO-20, XO-30, XO-40,
XO-50, XO-60, XO-77 and a
sterilization tray



X-Otome Kit K-7700

Includes one of each:
XO-77, XO-77AL, XO-77AR,
XO-77M, XO-77D and a
sterilization tray





D-LUX Elevators

ATRAUMATIC EXTRACTIONS
by A. Titan

Standard elevators are a thing of the past with the new D-Lux Elevators.



LUX 5

5mm



LUX 4

4mm




LUX 3

3mm




LUX 2

2mm



LUX 5A

5mm



LUX 4A

4mm



LUX 3A

3mm



LUX 2A


2mm



- A necessity to performing an atraumatic extraction
- Expedites patient healing time by preserving tissue and bone
- Sharp honed tips compress the alveolar & sever the PDL
- 100% USA 440A stainless steel
- D-Lux elevators are available individually in straight and curved patterns with millimeter widths of 2,3,4, and 5.

D-Lux Straight Kit K-5500S

Contains:
2mm, 3mm, 4mm, and 5mm
straight elevators with a
sterilization cassette



D-Lux Curved Kit K-5500A

Contains:
2mm, 3mm, 4mm, and 5mm
curved elevators with a
sterilization cassette



Complete Kit K-5500

Contains:
3 curved D-Lux elevators, 3
straight D-Lux elevators and
a sterilization cassette

