

# LEADER IN VETERINARY POINT OF CARE VISCOELASTIC COAGULATION MONITORING

OVEN C

## Why <mark>VCM Vet?</mark>

#### 1. Ease of Use

In just three steps VCM Vet™ will provide a patient's global hemostasis profile.

## 2. Flexibility

Whether running in the most up to date facilities or in the field, VCM Vet™ provides rapid, reproducible and high quality results.

#### 3. Affordable

More complete, actionable, real-time coagulation information in a cost-effective solution. VCM Vet's viscoelastic testing gives you a greater understanding of your patients' coagulation status.

#### 4. Patient Side Results

Immediately at your patient's side assess hemostasis status, diagnose and adjust therapies without waiting for results to come back from the lab. Real-time, point of care viscoelastic coagulation monitoring.

### 5. Reliablilty

Our automated, chemical activator-free test limits and prevents variation between operators, ensuring your results remain accurate regardless of user.

## Contact Us!

Email: info@vcmvet.com

@VCM\_Vet

Phone: 1-833-VCMVET





