

bili·swaddle™

neonatal procedural support garment

Little Sparrows
Technologies

Keep baby calm

4-way stretch fabric and patented peanut-shaped design encourages "hands over heart" position, keeping baby cocooned similar to being in the womb and reducing unnecessary waking from the startle reflex.

Provide developmental care

Ample leg room allows for sufficient flexion to promote proper hip development. A double zipper allows quick and easy access to patient diaper area and feet for pulse oximeter, ID band checks and heel stick blood sampling.



bili-swaddle™ was co-designed with our friends at Woombie® and is a Little Sparrows exclusive

Support treatment


The mesh fabric used in our swaddle allows 90% light transmittance to provide procedural support during phototherapy treatment.

Care providers also have easy access for IV site visualization and vital signs checks without overstimulating baby by wrapping and unwrapping swaddle blankets.

www.littlesparrowstech.com/phototherapy-bilihut

bili·swaddle™

when an extra hug
is needed



developmentally sensitive
procedural support during
phototherapy treatment