

# The Optimal Breast Center Workflow

## Take a walk through a breast center powered by MagView

MagView software streamlines and automates your mammography workflow for superior efficiency.



1

Patient receives appointment reminder text and completes history online prior to visit.

2

Patient completes history on tablet in waiting room, if needed.



3

Technologist verifies all information with patient and performs exam.



4

Radiologist reads and reports via structured reporting or 3rd party dictation or hybrid of both.



### Streamlined Reporting

Various integration options for the most efficient reading workflow:

- Luminary Workflow Manager
- AI Integrations
- PowerScribe Data Point Integration
- Pathology Interface
- Risk Assessment Module
- Workstation Sync

### ! Risk Assessment

- Facility chooses risk model
- Integration of automated breast density score available
- Patient's risk score flagged for technologist and radiologist
- Automatic recommendations for Ultrasound, MRI, or genetic testing

5



### Results

Patient receives lay results letter via online portal or through the mail.



### Reminders

Patient receives reminder texts before follow-up or annual appointment.

Patient returns for follow-up or annual appointment.



## ✓ Track Patients



Patient's needs tracked ongoing to ensure nothing is missed.

## ✓ Monitor Breast Center Performance



- Gain a comprehensive view
- Built-in reports and dashboards
- Adhere to regulatory requirements
- Prepare for MQSA inspection
- Enable industry accreditations